

Product Information

Justification:	Function:	Electronic signature justification:	Signed: Date (yyyy-mm-dd) - Time (hh:mm):
Issued:	QA	Adel Khwatmi - ADEKHW	2023-01-20 - 10:04
Reviewed:	QA	Abdallah Almashharawi - ABDALM	2023-01-20 - 10:55
Approved:	QA	Elin Andersson - ELIAND	2023-01-25 - 07:58
Released:	QA	Adel Khwatmi - ADEKHW	2023-02-07 - 08:52

This document has been electronically signed by the persons above.

TrachPhone®



Product description:

TrachPhone heats and humidifies the inhaled air and partially restores breathing resistance. It can be occluded with a finger to facilitate speech. After release the valve will open automatically. TrachPhone is connected to an ISO 15 tube. An integrated suction port makes it possible to clean the tracheostomy tube from mucus as needed. An oxygen tubing can be connected via the oxygen connector present on TrachPhone.

Product Information

Document ID: PF023-01-TechInfo **Edition:** 07

Manufacturer: Atos Medical AB
Kraftgatan 8
SE-242 35 Hörby, Sweden

Classification: (EU) (MDD 93/42/EEC) Class IIa (1.2 Rule 2)

Intended Use: For patients breathing spontaneously via an ET tube or a tracheostomy tube in the hospital or at home.

Use specifications: **Intended medical indication:** Product for rehabilitation for patients breathing through a tracheostoma.

Intended patient population: Patients of any age. Cognitive ability, by a clinician judged as sufficient. Manual dexterity, by a clinician judged as sufficient. Not intended for patients with mechanical ventilation. Not intended for patients with a low tidal volume.

Intended usage: Single use. Over the counter.

Intended part of the body/type of tissue applied to or interacted with:

The device will contact intact skin and mucosal membrane and as external communicating device the contact mode with tissue is indirect via air.

Intended user profile: The product is supposed to be handled by the patient but is also handled by physicians, trained nurses, SLPs, clinicians and caregivers.

Intended conditions of use:

Environment: Home use (normal daily environment without any or environmental restrictions regarding temperature, moisture etc.).
Outpatient clinic use. Hospital use.

Frequency of use: Continuous use.

Replacement rate: Max usage for 24 hours. Replacement is performed by the patient, clinician or caregiver.

Contraindications: Do not use beyond recommended tidal volume range, as the added dead space may cause CO₂ retention at too low tidal volumes. A too high tidal volume may lead to unsatisfactory humidification.
Do not use on dehydrated patients or patients with very heavy secretions from the lungs and airways.

CE Mark: Yes. Devices are CE-marked.

GMDN code: 58705 (Tracheostoma protective filter)

Sterilization: Non-Sterile

Raw material: Polypropylene (PP), thermoplastic elastomers (TPE) and polyurethane (PUR).

Latex information: Not manufactured with natural rubber latex.

Biological origin: The device is not manufactured with materials derived from human or animal source.

Handling and storage: Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.

Product Information

Waste handling and disposal:

Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.

Hazardous components:

None

Expiration date:

3 years after manufacturing.

Packaging:

TrachPhone is available as 50, 30 and 5 pack.
 Each TrachPhone is packed in a plastic bag.
 50 / 30 / 5 plastic bags are packed in an inner box (a total of 50 /30 / 5 cassettes).

Devices under Basic UDI-DI: 7331791-HME-0-000-0006-XT

REF	Name	UDI-DI
7704	TrachPhone (50 pcs)	07331791002861
7707	TrachPhone (30 pcs)	07331791009693
7723	TrachPhone (5 pcs)	07331791015854

Atos Medical AB compatible products:

Range	BASIC UDI-DI
Fits on standard 15 mm ISO connector.	N/A

Freevent® DualCare



Product description:

Freevent DualCare Speaking valve is a valve with a silicone membrane and a rotatable lid. HME DigiTop is a top that can be occluded with two digits to enable speech. Both these speaking devices are attached to either Freevent HME 15 Regular or Freevent HME 22 Regular before use.

Freevent Connection Strap is a clip with a string that is used to secure Freevent DualCare to the patient's neckband.

Removal aid is a plastic clamp that is pressed together by finger force to clamp the HME at HME removal from the speaking devices.

Document ID: PF068-01-TechInfo

Edition: 2.0

Manufacturer: Atos Medical AB
Kraftgatan 8
SE-242 35 Hörby, Sweden

Classification: (EU) Class I (Rule 1)
2017/745

Intended Use: **Freevent DualCare** is a combined Speaking Valve and Heat and Moisture Exchanger (HME) intended for spontaneously breathing tracheostomized patients using a tracheostomy tube with a deflated cuff, or a tracheostomy tube without cuff.

In HME-mode the device conditions inhaled air by retaining heat and moisture from the exhaled air. By turning the lid of the Speaking Valve into speaking mode air is re-directed to enable speech.

The entire device is for single patient use and the HME-part is for single use.

Freevent HME 15 Regular is a Heat and Moisture Exchanger (HME) intended for spontaneously breathing tracheostomized patients using a tracheostomy tube with an ISO 15 mm connector. The HME conditions inhaled air by retaining heat and moisture from the exhaled air. The device also partially restores lost breathing resistance.

The HME is used in combination with Freevent DualCare Speaking valve/Freevent DualCare Speaking Valve Blue, with Freevent HME DigiTop/Freevent HME DigiTop Blue, or with HME DigiTop O2.

The HME is for single use, i.e. it has to be exchanged at least every 24 hours.

Freevent HME 22 Regular is a Heat and Moisture Exchanger (HME) intended for spontaneously breathing tracheostomized patients using a tracheostomy tube with a Ø22mm connector. The HME conditions inhaled air by retaining heat and moisture from the exhaled air. The device also partially restores lost breathing resistance.

The HME is used in combination with Freevent DualCare Speaking Valve/Freevent DualCare Speaking Valve Blue, with Freevent HME DigiTop/Freevent HME DigiTop Blue, or with HME DigiTop O2. The HME is for single use, i.e. it has to be exchanged at least every 24 hours.

Use specifications: Freevent DualCare

Intended medical indication

Product for rehabilitation for patients breathing through a tracheostoma.

Intended patient population

Patients of any age.

Cognitive ability, by a clinician judged as sufficient.

Manual dexterity, by a clinician judged as sufficient.

Not intended for patients with mechanical ventilation.

Not intended for patients with a low tidal volume.

Intended usage

HME: Single use, Prescription only .

Speaking Valve: Single patient multiple use, Prescription only.

Removal Aid and Connection Strap: Single patient multiple use, Over-the-counter.

Intended part of the body/type of tissue applied to or interacted with

The device will contact intact skin and mucosal membrane and as external communicating device the contact mode with tissue is indirect via air.

Intended user profile

The product is supposed to be handled by the patient but is also handled by physicians, trained nurses, SLPs, clinicians and caregivers.

Intended conditions of use

Environment: Home use (normal daily environment without any or environmental restrictions regarding temperature, moisture etc.).

Outpatient clinic use. Hospital use.

Frequency of use: Continuous use.

Replacement rate: Max usage for 24 hours for HME. Replacement is performed by the patient, clinician or caregiver.

Product Information

Freevent HME 15 Regular and Freevent HME 22 Regular

Intended medical indication

Product for rehabilitation for patients breathing through a tracheostoma.

Intended patient population

Patients of any age.

Cognitive ability, by a clinician judged as sufficient.

Manual dexterity, by a clinician judged as sufficient.

Not intended for patients with mechanical ventilation.

Not intended for patients with a low tidal volume.

Intended usage

Single use, Prescription only.

Intended part of the body/type of tissue applied to or interacted with

The device will contact intact skin and mucosal membrane and as external communicating device the contact mode with tissue is indirect via air.

Intended user profile

The product is supposed to be handled by the patient but is also handled by physicians, trained nurses, SLPs, clinicians and caregivers.

Intended conditions of use

Environment: Home use (normal daily environment without any or environmental restrictions regarding temperature, moisture etc.).

Outpatient clinic use. Hospital use.

Frequency of use: Continuous use.

Replacement rate: Max usage for 24 hours. Replacement is performed by the patient, clinician or caregiver.

Contraindications: Freevent DualCare

General

Freevent HME 15 or 22 in combination with either Freevent DualCare Speaking Valve or Freevent HME DigiTop are contraindicated for: Use in combination with an in-line ventilator.

- Patients without the physical, cognitive, or mental ability required to attach, remove, or operate the devices themselves, should not use the devices independently and should only use them if they are under sufficient supervision of a clinician or a trained caregiver.
- The devices should not be used by patients with a low tidal volume, as the added dead space may cause CO₂ (Carbon dioxide) retention.
- Unresponsive or sedated patients.

The HME 15 or 22 in combination with either Speaking Valve or HME DigiTop must NOT be used on a single lumen tube (tube without an inner tube), unless the patient or caregiver is able to reinsert the tube themselves after accidental dislodgment or emergency replacement.

Speaking valve specific

The use of Speaking Valve (in combination with HME 15 or 22) is additionally contraindicated for the following patient groups:

- Laryngectomized patients since the device will prevent the ability to exhale if the Speaking Valve is unintentionally set to speaking mode.
- Patients suffering from severe aspiration.
- Patients with severe laryngeal or upper airway obstruction such as significant tracheal and/or laryngeal stenosis, since this may cause air trapping.
- Patients with very thick and copious secretions which might block the device.

DO NOT use the Speaking Valve:

- In combination with a tracheostomy tube with the cuff inflated. The cuff must be completely deflated before placing and during all use of the Speaking Valve.
- In combination with a tracheostomy tube with a foam cuff.
- In combination with a tracheostomy tube with a self-inflating cuff.
- When the size of the tracheostomy tube does not allow for airflow through the upper airways.
- In combination with an endotracheal tube.

Use of the Speaking Valve in these circumstances can restrict exhalation through the upper airways and cause suffocation!

DO NOT use the Speaking Valve during sleep since the airway could be blocked unintentionally. During sleep the HME DigiTop (in combination with HME 15 or 22) should be used instead.

Product Information

Freevent HME 15 Regular and Freevent HME 22 Regular

Freevent HME 22 in combination with either Freevent DualCare Speaking Valve or Freevent HME DigiTop is contraindicated for:

- Use in combination with an in-line ventilator.
- Patients who are unable to handle or remove the device themselves when needed, unless the patient is under constant supervision of a clinician or a trained caregiver. For example, patients who are unable to move their arms, patients with decreased levels of consciousness, or patients with diseases that put them at risk for unpredictable periodic loss of consciousness.
- Patients who cannot tolerate the added dead space of 5 ml, or who cannot tolerate the added breathing resistance of 170 Pa / 1.7 cm H₂O. This should be evaluated by a clinician.
- Unresponsive or sedated patients. The patient must be responsive and attempting to communicate in order to use the device. The patient should be able to follow instructions.

The HME 22 in combination with either Speaking Valve or HME DigiTop must NOT be used on a single lumen tube (tube without an inner tube), unless the patient or caregiver is able to reinsert the tube themselves after emergency removal.

CE Mark:	Yes. Devices are CE-marked.
GMDN code:	58705 (Tracheostoma protective filter) for REF 7742, 7745 and 7747. 36071 (Tracheostomy tube speech valve) for REF 7740, 7741, 7744, 7746 and 7755.
Sterilization:	Non-Sterile
Raw material:	Freevent Speaking Valve: PP, Silicone, and POM. Freevent Speaking Valve Blue: PP, Silicone, and POM. Freevent DigiTop: POM. Freevent HME 15 Regular: POM, HDPE, and Polyester-based Polyurethane foam. Freevent HME 22 Regular: Polypropylene (PP) with white masterbatch, and Polyester-based Polyurethane foam. Freevent Connection Strap: Polyester braided suture, POM, and PP. Removal Aid: POM.
Latex information:	Not manufactured with natural rubber latex.
Biological origin:	The device is not manufactured with materials derived from human or animal source.
Handling and storage:	Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.
Waste handling and disposal:	Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.
Hazardous components:	None.
Expiration date:	3 years after manufacturing.

Product Information

Packaging:

Product	Contents
7740 Freevent DualCare Set 22	HME 22 Regular, 30 pcs (10 pcs/bag) are packed in plastic bag of polyethylene. Speaking Valve, 1 pc is packed in a plastic jar of polypropylene. HME DigiTop, 1 pc is packed in a plastic bag of polyethylene. Removal Aid, 1 pc is packed in a plastic bag polyethylene. Connection Strap, 1 pc is packed in a plastic bag of polyethylene. The products, 2 pcs warning label sheets and instructions for use are packed in a cardboard box.
7741 Freevent DualCare Set 15	HME 15 Regular, 30 pcs (10 pcs/bag) are packed in plastic bag of polyethylene. Speaking Valve, 1 pc is packed in a plastic jar of polypropylene. HME DigiTop, 1 pc is packed in a plastic bag polyethylene. Removal Aid, 1 pc is packed in a plastic bag polyethylene. Connection Strap, 1 pc is packed in a plastic bag polyethylene. The products, 2 pcs warning label sheets and instructions for use are packed in a cardboard box.
7742 Freevent HME 15 Regular (30pcs)	HME 15 Regular, 30 pcs (10 pcs/bag) are packed in plastic bag of polyethylene. The products are packed in a cardboard box.
7744 Freevent DualCare Speaking Valve	Speaking Valve, 1 pc is packed in a plastic jar of polypropylene. HME DigiTop, 1 pc is packed in a plastic bag polyethylene. Connection Strap, 1 pc is packed in a plastic bag polyethylene. The products, 2 pcs warning label sheets and instructions for use are packed in a cardboard box.
7745 Removal Aid	Removal Aid, 1 pc is packed in a plastic bag polyethylene. The product and instructions for use are packed in a bubble plastic bag.
7746 Freevent Connection strap	Connection Strap, 2 pcs (1 pc/bag) are packed in plastic bag polyethylene. The products and instructions for use are packed in a bubble plastic bag
7747 Freevent HME 22 Regular (30pcs)	HME 22 Regular, 30 pcs (10 pcs/bag) are packed in plastic bag of polyethylene. The products are packed in a cardboard box.
7755 Freevent DualCare Speaking Valve Blue	Speaking Valve Blue, 1 pc is packed in a plastic jar of polypropylene. HME DigiTop Blue, 1 pc is packed in a plastic bag polyethylene. Connection Strap, 1 pc is packed in a plastic bag polyethylene. The products, 2 pcs warning label sheets and instructions for use are packed in a cardboard box.

Devices under Basic UDI-DI: 7331791-HME-0-000-0005-XQ

REF	Name	UDI-DI
7740	Freevent DualCare Set 22	7331791015038
7741	Freevent DualCare Set 15	7331791015021
7744	Freevent DualCare Speaking Valve	7331791015045
7745	Removal Aid	7331791008221
7746	Freevent Connection strap	7331791008238
7755	Freevent DualCare Speaking Valve Blue	7331791015052

Devices under Basic UDI-DI: 7331791-HME-0-000-0010-XE

REF	Name	UDI-DI
7742	Freevent HME 15 Regular (30pcs)	7331791015069

Devices under Basic UDI-DI: 7331791-HME-0-000-0011-XH

REF	Name	UDI-DI
7747	Freevent HME 22 Regular (30pcs)	7331791015076



HME 15



HME 22



Speaking Valve (Blue)



DigiTop (Blue)



Connection strap



Removal aid

Atos Medical AB compatible products:

Range	BASIC UDI-DI
HME DigiTop O2	7331791-HME-A-000-0007-FH

Document Approvals
Approved Date: 2024-05-23

Task: Approval Task Verdict: Approve	ABDALM Abdallah Almashharawi, Sustaining Engineer (abdallah.almashharawi-atosmedical@coloplast.com) Issuer 21-May-2024 07:22:27 GMT+0000
Task: Approval Task Verdict: Approve	SOFTHO Sofia Thomasson, Regulatory Affairs Professional (sofia.thomasson-atosmedical@coloplast.com) Quality 21-May-2024 07:44:02 GMT+0000
Task: Final Approval Verdict: Approve	KARNIL Karolina Nilsson, Head of Regulatory Affairs (karolina.nilsson-atosmedical@coloplast.com) Regulatory 23-May-2024 12:49:44 GMT+0000

Technical Info / Material Data Sheet

Document ID:
 PF068-07-Tech Info

Edition: 00

REF Number 7756

Product Name HME DigiTop O2 (REF7756)

Models: One variant, fitting for 22mm HME Cassette.
 One product, each containing one DigiTop O2 + Instructions For Use.

Classification: Class IIa, 1.2 rule 2
 (MDD 93/42/EEC)

CE Mark: Yes

GMDN code: 58705

Produced by: Atos Medical AB
 Kraftgatan 8
 P.O. Box 183
 242 22 Hörby
 Sweden

Intended Use: The HME DigiTop O2 is an accessory to ProTrach HMEs and Provox FreeHands HME's. For patients spontaneously breathing through a tracheostoma and having a need of extra oxygen.

Description: HME DigiTop O2 is a top that can be occluded with two digits to enable speech. The device shall be attached to either ProTrach HME 15 or HME 22 before use. The oxygen connector port on the device shall be connected to an oxygen supply via a tube

Sterilization: Non-sterile

Raw material: HME DigiTop O2: Blue POM .

Latex information The device is not manufactured with natural rubber latex.

Biological origin: The device is not manufactured with any materials derived from human or animal source.

Handling and storage: Keep dry and away from sunlight. Temperature limit: 2-42 °C.

Waste handling and disposal: Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.

Hazardous components: None.

Expiration date: 3 years after manufacturing.

Packaging: REF 7756, HME DigiTop O2: – box with 1 pcs plastic jar with 1 pc HME DigiTop O2 + 1 pc IFU REF 10721.

Technical Info / Material Data Sheet

Reviewed by:


.....
Vice President QA&RA


.....
Date

Approved by:


.....
Vice President Design Control


.....
Date

Justification:	Function:	Electronic signature justification:	Signed: Date (yyyy-mm-dd) - Time (hh:mm):
Issued:	QA	Carolina Johansson - SEHRBJNC	2022-04-21 - 07:07
Reviewed:	QA	Karolina Nilsson - KARNIL	2022-04-21 - 09:54
Approved:	DD	Diana Tieger - DIATIE	2022-04-21 - 10:41
Released:	QA	Carolina Johansson - SEHRBJNC	2022-05-19 - 15:05

This document has been electronically signed by the persons above.

Freevent® XtraCare and Freevent® XtraCare Mini



Figure 1- FreeVent XtraCare Blue



Figure 2- FreeVent XtraCare Mini Blue

Product description:

Freevent XtraCare and Freevent XtraCare Mini are Heat and Moisture Exchangers combined with an electrostatic filter (HMEF). The HME is impregnated with a hygroscopic salt and conditions the inhaled air. The electrostatic filter reduces the inhalation of particles such as viruses, bacteria, pollen and other particulate matter through the tracheostoma. Freevent XtraCare and Freevent XtraCare Mini have a 15 mm ISO connector for connection to a tracheostomy tube.

Freevent XtraCare comes in two colors, white and blue. Each color comes in two package sizes, 5 pcs and 30 pcs.

Freevent XtraCare Mini comes in three colors, white, blue and pink. Each color comes in a package size of 30 pcs. Freevent XtraCare Mini White come in an additional package size of 5 pcs.

Freevent XtraCare and XtraCare Mini can be connected to oxygen tubing using the Freevent O2 Adaptor respectively O2 Adaptor mini (accessory).

Product Information

Document ID:	PF069-01-TechInfo	Edition:	05
Manufacturer:	Atos Medical AB Kraftgatan 8 SE-242 35 Hörby,Sweden		
Classification: (EU) 2017/745	Class I (Rule 1)		
Intended Use:	Freevent® XtraCare and Freevent® XtraCare Mini are single use Heat and Moisture Exchangers with electrostatic filters (HMEF) that condition and filter inhaled air in patients spontaneously breathing through a tracheostoma		
Use specifications:	<p>Intended medical indication: Patients breathing through a tracheostoma, long-term and short-term, independent of underlying condition. Especially intended for patients with a need for enhanced protection against microorganisms/pathogens, pollen, and other particles.</p> <p>Intended patient population: For patients with any health condition who breathe spontaneously through a tracheostoma and can tolerate the added dead space of the product and the added breathing resistance.</p> <p>Intended usage: Disposable single use product. Can be used 24/7 and shall be replaced if the breathing resistance has become too high e.g. when saturated with mucus, or if the 24 hours limit has been reached. Prescription only.</p> <p>Intended part of the body/type of tissue applied to or interacted with: To be applied on a tracheostomy tube or similar device with a 15mm connector.</p> <p>Intended user profile: Clinicians, caregivers, patient, depending on the condition of the patient.</p> <p>Intended conditions of use: Environment of use: Hospitals, ICU, Sub-acute care institutions, and home, indoors and outdoors. It does not affect the patient's mobility.</p>		
Contraindications:	<p>patients who:</p> <ul style="list-style-type: none"> • are under any form of mechanical ventilation. • are unable to handle or remove the device themselves when needed, and who are not under constant supervision of a clinician or a trained caregiver. • cannot tolerate the added dead space. 		
CE Mark:	Yes. Devices are CE-marked.		
GMDN code:	58705 (Tracheostoma protective filter)		
Sterilization:	Non-Sterile		
Raw material:	<p>Plastic parts (Base and Housing): Polypropylene (PP) with white, blue or pink PP masterbatch.</p> <p>Foam: Polyurethane (PUR) with Calcium Chloride (CaCl₂)</p> <p>Filter: Acrylic fiber attached to Polypropylene (PP) spunbonded scrim</p>		
Latex information:	Not manufactured with natural rubber latex		

Product Information

Biological origin:	The device is not manufactured with materials derived from human or animal source.
Handling and storage:	Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.
Waste handling and disposal:	Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.
Hazardous components:	None
Expiration date:	3 years after manufacturing.
Packaging:	<p>Freevent XtraCare 5pcs/30pcs (1pc/bag) are packed in plastic bag of polyethylene. The products and instructions for use are packed in a cardboard box.</p> <p>Freevent XtraCare Mini 5pcs/30pcs (1pc/bag) are packed in plastic bag of polyethylene. The products and instructions for use are packed in a cardboard box.</p>

Devices under Basic UDI-DI: 7331791-HME-0-000-0004-XM

REF	Name	UDI-DI
7767	Freevent XtraCare, white (30 pcs)	07331791008948
7768	Freevent XtraCare, blue (30 pcs)	07331791008955
7789	Freevent XtraCare, white (5 pcs)	07331791008962
7788	Freevent XtraCare, blue (5 pcs)	07331791008979
8004	Freevent XtraCare Mini white (30 pcs)	07331791014901
8005	Freevent XtraCare Mini blue (30 pcs)	07331791014918
8006	Freevent XtraCare Mini pink (30 pcs)	07331791014925
8008	Freevent XtraCare Mini white (5 pcs)	07331791014932

Atos Medical AB compatible products:

Range	BASIC UDI-DI
Freevent O2 Adaptor	7331791-HME-A-000-0001-EX
Freevent O2 Adaptor mini	7331791-HME-A-000-0001-EX
Provox BasePlate Adaptor	7331791-HME-A-000-0003-F5

Justification:	Function:	Electronic signature justification:	Signed: Date (yyyy-mm-dd) - Time (hh:mm):
Issued:	QA	Adel Khwatmi - ADEKHW	2022-11-25 - 15:02
Reviewed:	QA	Abdallah Almashharawi - ABDALM	2022-11-25 - 15:23
Approved:	QA	Elin Andersson - ELIAND	2022-11-28 - 11:27
Released:	QA	Abdallah Almashharawi - ABDALM	2022-12-14 - 10:10

This document has been electronically signed by the persons above.

Product Information

Freevent® O2 Adaptors



Product description:

Freevent O2 Adaptors are accessories that fit Freevent XtraCare and Freevent XtraCare Mini. They are clicked over the base of the HME and the combined device is attached to the patient’s tracheostomy tube, or similar device. Additional oxygen can then be supplied via the oxygen port of the O2 Adaptor. Freevent O2 Adaptors are single use devices and should be replaced if they become dirty, or at least every 24 hours.



The O2 Adaptor mounted on a Blue Freevent XtraCare

Document No: 10000038367 Edition: 04 Release date: 2022-12-14
Released

Product Information

Document ID:	PF069-02-TechInfo	Edition:	04
Manufacturer:	Atos Medical AB Kraftgatan 8 SE-242 35 Hörby,Sweden		
Classification: (EU) 2017/745	Class IIa (1.2 rule 2)		
Intended Use:	Freevent O2 Adaptor and Freevent O2 Adaptor Mini are single use accessories used together with Freevent XtraCare or Freevent XtraCare Mini respectively. The devices are intended to enable additional oxygen supply for patients breathing through a tracheostoma during use of Freevent XtraCare and Freevent XtraCare Mini.		
Use specifications:	<p>Intended medical indication Patients breathing through a tracheostoma, long-term and short-term, independent of underlying condition.</p> <p>Intended patient population For patients with any health condition who breathe spontaneously through a tracheostoma and use applicable Freevent product.</p> <p>Intended usage Disposable single use product. Should be changed at least every 24 hours. For prescription only.</p> <p>Intended part of the body/type of tissue applied to or interacted with: To be applied on applicable Freevent XtraCare, which in turn is connected to a tracheostomy tube or similar device.</p> <p>Intended user profile Clinicians, caregivers, patients, depending on the condition of the patient.</p> <p>Intended conditions of use Environment of use: Hospitals, ICU, Sub-acute care institutions, and home, indoors and outdoors. The device itself does not affect the patient's mobility.</p>		
Contraindications:	No known contraindications.		
CE Mark:	Yes. Devices are CE-marked.		
GMDN code:	58705 (Tracheostoma protective filter)		
Sterilization:	Non-Sterile		
Raw material:	Polypropylene (PP)		
Latex information:	Not manufactured with natural rubber latex		
Biological origin:	The device is not manufactured with materials derived from human or animal source.		
Handling and storage:	Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.		

Product Information

Waste handling and disposal:	Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.
Hazardous components:	None
Expiration date:	3 years after manufacturing.
Packaging:	Freevent O2 Adaptors are single packed in a plastic bag of polyethylene and then 10 pieces in a cardboard box together with IFU.

Devices under Basic UDI-DI:

REF	Name	UDI-DI
7769	Freevent O2 Adaptor 10pcs	07331791008986
8007	Freevent O2 Adaptor Mini 10pcs	07331791015311

Atos Medical AB compatible products:

Range	BASIC UDI-DI
Freevent XtraCare White/Blue	7331791-HME-0-000-0004-XM
Freevent XtraCare Mini White/Blue/Pink <i>(for O2 adaptor Mini)</i>	7331791-HME-0-000-0004-XM