

# PROVOX Life

## LaryButton and HMEs

CE 0413 LaryButton

CE HME

UK  
CA HME

LaryButton

**Rx**  
**ONLY**

LaryButton

HME


2 °C  
36 °F


42 °C  
108 °F

RT


## Provox Life LaryButton


1


2


3


3 a-1


3 a-2


3 b-1


3 b-2

3 b-3


3 b-4


4a


4b


4c


5a


5b


5c


5d

5e


5f


5g


Provox Life Home HME/Go HME/  
Energy HME/Protect HME

6.


7.

8.


## Provox Life Night HME

9.


10.





11.

12.


13.

17.

18.

## Prescription information

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## Disclaimer

Atos Medical offers no warranty—neither expressed nor implied—to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## Contents

EN - ENGLISH LARYBUTTON.....	16
EN - ENGLISH HME .....	25
DE - DEUTSCH LARYBUTTON .....	30
DE - DEUTSCH HME .....	40
NL - NEDERLANDS LARYBUTTON .....	46
NL - NEDERLANDS HME .....	56
FR - FRANÇAIS LARYBUTTON.....	62
FR - FRANÇAIS HME.....	71
IT - ITALIANO LARYBUTTON.....	77
IT - ITALIANO HME.....	86
ES - ESPAÑOL LARYBUTTON .....	91
ES - ESPAÑOL HME.....	100
PT - PORTUGUÊS LARYBUTTON.....	105
PT - PORTUGUÊS HME.....	114
SV - SVENSKA LARYBUTTON.....	119
SV - SVENSKA HME.....	128
DA - DANSK LARYBUTTON .....	133
DA - DANSK HME .....	142
NO - NORSK LARYBUTTON .....	147
NO - NORSK HME .....	156

FI - SUOMI LARYBUTTON.....	161
FI - SUOMI HME.....	170
IS - ÍSLENSKA LARYBUTTON.....	175
IS - ÍSLENSKA HME.....	184
ET - EESTI LARYBUTTON.....	189
ET - EESTI HME.....	198
CS - ČESKY LARYBUTTON .....	203
CS - ČESKY HME.....	212
HU - MAGYAR LARYBUTTON .....	218
HU - MAGYAR HME .....	227
PL - POLSKI LARYBUTTON.....	232
PL - POLSKI HME.....	242
HR - HRVATSKI LARYBUTTON .....	261
HR - HRVATSKI HME .....	256
EL - ΕΛΛΗΝΙΚΑ LARYBUTTON .....	261
EL - ΕΛΛΗΝΙΚΑ HME.....	272
JA - 日本語 LARYBUTTON.....	278
JA - 日本語 HME.....	287
AR - العربية LARYBUTTON.....	298
AR - العربية HME.....	299

## **Provox Life System**

### **1. Description of the device**

Provox Life devices are designed to work together as a system. Provox Life System consists of Provox Life Adhesives, Provox Life HMEs, Provox Life LaryTubes, Provox Life LaryButtons and Provox Life accessories range of devices.

### **2. CONTRAINDICATIONS**

The devices should only be used in accordance with the Instructions for Use.

Patients without the physical, cognitive, or mental ability required to attach, remove or operate the devices themselves, should not use the devices independently and should only use them if they are under sufficient supervision of a clinician or a trained caregiver.

#### Only valid for HMEs

The devices should not be used by patients with a low tidal volume, as the added dead space may cause CO<sub>2</sub> (Carbon dioxide) retention, see 2.1. Technical data.

## **Provox Life LaryButton**

### **1. Intended use**

Provox Life LaryButton is a single patient use self-retaining holder for Provox Life HMEs and accessories after total laryngectomy.

For laryngectomized patients with a shrinking

tracheostoma it is also used to maintain the tracheostoma for breathing.

## 2. CONTRAINDICATIONS

Provox Life LaryButton is not intended to be used by patients that:

- are under any form of mechanical ventilation.
- have damaged tracheostoma tissue.

## 3. Description of the device

Provox Life LaryButton is a self-retaining tracheostoma button made of medical grade silicone rubber. The device is delivered single packed, non-sterile and ready to use. The purpose of the device is to create a comfortable and airtight fit between Provox Life LaryButton and the tracheostoma, and also to provide attachment for Provox Life HMEs and Provox Life accessories.

The different parts of Provox Life LaryButton are (Fig. 1):

- a) Retention Collar
- b) Shaft
- c) Shield (conical)
- d) HME and Accessory Holder
- e) Wings

## 4. WARNINGS

- Always use Provox Life LaryButton of an appropriate size. If Provox Life LaryButton is too wide or too long it may cause tissue damage or irritation. If it is too long it may also cause dislodgement of the voice prosthesis or Provox XtraFlange (Fig. 3 a-2). Furthermore,

insertion or removal of the HME as well as occluding the stoma by pressing the lid of the HME to speak may exert pressure on the voice prosthesis if Provox Life LaryButton is too long.

- Always make sure that the shield (Fig. 1, c) is large enough in relation to the size of the patient's stoma. If the product is pushed too far into the trachea it may obstruct breathing.
- DO NOT cut off the wings or otherwise mechanically modify (e.g. fenestrate) Provox Life LaryButton since it will change the outer diameter and/or the mechanical stability of the product. This may result in Provox Life LaryButton being pushed into the trachea, which may obstruct breathing. A modified Provox Life LaryButton may also cause irritation and/or bleeding.
- DO NOT reuse Provox Life LaryButton between patients. The device is intended for single patient use only. Reuse between patients may cause cross contamination.
- DO NOT use force and do not fold the part of Provox Life LaryButton that holds the HME during insertion. If the product is pushed too far into trachea it may obstruct breathing and the use of force during insertion may cause bleeding.
- DO NOT push the tip of the hemostats or forceps (if used) through Provox Life LaryButton (Fig. 3 b-2). If the instrument tip is protruding from the button during insertion, the tissue of the trachea could be harmed and/or the voice prosthesis could become dislodged and subsequently be aspirated or ingested.

- DO NOT use lubricating gels if the patient has allergies related to these substances.
- DO NOT at any time attempt to clean the device while inside the stoma since this could cause tissue damage and dislodgement of voice prosthesis.
- Make sure that the patient has been trained in the use of the device, and has demonstrated the ability to understand and consistently follow the Instructions for Use without clinician supervision.

## 5. PRECAUTIONS

- DO NOT use a dirty or contaminated Provox Life LaryButton. Clean and disinfect the device according to the cleaning and disinfection instructions below.
- DO NOT use toothed or sharp instruments at any time during insertion since they could damage Provox Life LaryButton and cause tissue damage.
- DO NOT use a damaged device since this might cause aspiration of small parts or cause formation of granulation tissue.
- DO NOT use oil-based lubricants (e.g., Vaseline), they may weaken, damage, or destroy the product. Only use water-soluble lubricants.
- Carefully monitor the tracheostoma tissue during radiation therapy and stop using the device if the stoma gets irritated or starts to bleed.
- Patients with bleeding disorders or patients on anticoagulants should not use the device if it causes re-occurring bleeding.
- Stop using the device if granulation tissue forms in the tracheostoma.

## 6. Instructions for use

### 6.1 Preparation

#### **Choose correct size**

- Choose the correct diameter and length of Provox Life LaryButton.
- To achieve an airtight fit and good retention the diameter of the Retention Collar (Fig. 1, a) of Provox Life LaryButton should exceed the diameter of the tracheostoma entrance.
- Ensure that the Retention Collar of Provox Life LaryButton fits right behind the edge of the tracheostoma. However if the device is too short it may fall out of the tracheostoma.

**CAUTION:** During fitting always make sure that Provox Life LaryButton does not reach the tracheal flange of the voice prosthesis. Optimal fit and performance of the Provox Life LaryButton is achieved if the patient has a so-called “tracheostoma lip” that is situated uniformly around the tracheostoma.

#### **Prepare for insertion**

- Always inspect the product before each use. Do not use if damaged or dirty.
- Make sure that your hands are clean and dry.
- If needed, lightly lubricate the tube with a small amount of the following recommended water-soluble lubricants: KY Jelly, SurgiLube and Xylocaine (lidocaine).

## 6.2 Operating instruction

### Insertion (Fig. 3):

#### **Alternative 3a:**

1. Fold the Retention Collar with your fingers (Fig. 3) and hold it firmly folded.
2. Insert the Retention Collar into the stoma and release it (Fig. 3 a-1 and 3 a-2).

#### **Alternative 3b:**

1. Fold the Retention Collar with your fingers (Fig. 3).
2. Use a pair of blunt forceps or a lockable blunt hemostat from the inside to keep it folded (Fig. 3 b-1 and b-2).  
**CAUTION: DO NOT** push the tip of the hemostats or forceps through Provox Life LaryButton (Fig. 3 b-2). If the instrument tip is protruding from the button during insertion, the tissue of the trachea could be harmed and/or the voice prosthesis could become dislodged and subsequently be aspirated or ingested.
3. Insert the Retention Collar into the stoma (Fig. 3 b-3).
4. Release the forceps to unfold the button (Fig. 3 b-4).

Check the correct fit by slightly pulling and/or rotating the button. Gently insert Provox Life HME/accessory into the HME and Accessory Holder.

**NOTE:** Some patients might experience more coughing after inserting Provox Life LaryButton and when touching the device to speak. This usually diminishes when the user gets used to the device. If discomfort occurs Provox Life LaryButton should be removed.

### **Device removal**

1. Hold Provox Life LaryButton in place and gently remove Provox Life HME/accessory from the HME and Accessory Holder (Fig. 4a).
2. Hold Provox Life LaryButton and gently turn it (about half a turn) to allow the Retention Collar to compress thereby reducing its width (Fig. 4b).
3. Remove the device by gently tilting it out of the stoma (Fig. 4c). Alternatively, the device can be removed by pushing with a finger on the shaft from one direction and tilt the device out.

**NOTE:** The insertion or removal procedures may occasionally cause some slight bleeding, irritation or coughing. Inform your clinician if the bleeding persists.

## **6.3 Cleaning and disinfection**

Clean and disinfect Provox Life LaryButton before each use. Only use drinking water to clean and rinse the device(Fig. 5a-5g).

1. Rinse Provox Life LaryButton with water (Fig. 5a).
2. Place Provox Life LaryButton in water-based mild soap for 15 minutes (Fig. 5b).
3. Clean inside of Provox Life LaryButton with Provox TubeBrush or Provox Swab and rinse with water (Fig. 5c + 5d).
4. Disinfect Provox Life LaryButton at least once a day with one of the following methods: Ethanol 70% for 10 minutes or Isopropylalcohol 70% for 10 minutes or Hydrogenperoxide 3% for 60 minutes (Fig. 5e).

5. Let Provox Life LaryButton dry for at least 2 hours or until it looks dry (Fig. 5f).
6. Rinse Provox TubeBrush or Provox Swab in water after use (Fig. 5g).

## 7. Storage instructions

Store Provox Life LaryButton in a clean and dry container at room temperature. Protect from direct sunlight.

## 8. Device lifetime and Disposal

Provox Life LaryButton may be used for a maximum period of 6 month. If the device shows any signs of damage it shall be replaced earlier.

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## 9. Compatible Products

Use only original components. Other products may cause product damage and/or malfunction.

Devices that can be used with Provox Life LaryButton are:

- Provox Life HMEs (Heat and Moisture Exchangers) for patients breathing through a tracheostoma.
- Provox Life Shower is used to prevent water from entering Provox Life LaryButton during showering.
- Provox Life BasePlate Adaptor, to attach products with an ISO 15 mm standard connector to Provox Life LaryButton.
- Provox TubeHolder or Provox LaryClip can be used to keep Provox Life LaryButton in place.

## 10. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

# Provox Life HME

## 1. Intended use

Provox Life HMEs are single-use heat and moisture exchangers for patients breathing through a tracheostoma.

## 2. Description of the device

Provox Life HMEs are single-use devices for pulmonary rehabilitation. They come in different levels of humidification, breathing resistance and filtration that makes them suitable for different situations.

The different Provox Life HMEs are:

- Home HME: when taking it easy
- Go HME: when you are out and about
- Energy HME: when physically active
- Protect HME: when you need protection from bacteria, virus, dust and pollen
- Night HME: when sleeping.

### 2.1 Technical data

	Pressure Drop at 30 l/min*(Pa)	Pressure Drop at 60 l/min*(Pa)	Pressure Drop at 90 l/min*(Pa)	Moisture loss at VT=1000 ml* (mg/l)
Home HME	60	200	450	19.5
Go HME	30	100	230	22.5

Energy HME	15	50	110	23
Protect** HME	55	180	320	23
Night HME	65	210	470	18.5

HME internal volume (Dead space)\* max 13 ml

\* According to ISO 9360

\*\* Bacterial Filtration Efficiency >98%, Virus Filtration Efficiency >98%. Test method adapted from ASTM F2101.

### 3. Warnings

- Do not exert unintentional physical pressure on the device. Unintentional or accidental closing or compression of the device may obstruct breathing.
- Do not squeeze, or apply excessive force during attachment or when the device is in place. Excessive force can result in the device being pushed through the adhesive coupling, and may obstruct breathing and block airways.
- Do not use non-compatible devices or non-original devices as it may cause personal injury or damage to the device. Do not dismantle or alter the device.
- Do not re-attach a device clogged with mucus, as it may lead to obstructed breathing.
- Do not attach a device that appears slippery by mucus or other fluids, as it will be hard to attach in a correct position.

### 4. Precautions

- Always test the function of the device prior to use. The lid should immediately return to its open position after releasing it. (Not applicable for Provox Life Night HME).

- Pay attention when the device needs to be removed or attached during the night as the night-time situation could make it more difficult to orient the device in a correct position, and it may also be difficult to detect if contaminated.
- To reduce the risk of water entering the airways, use a Provox Life Shower when showering.
- The device shall not be used beyond 24 hours, is not reusable and shall not be rinsed in water or in any other solutions. This will substantially reduce the function of the device and can increase the risk for infections due to bacterial colonization.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device. The device will become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.

## 5. Instructions for use

For guidance see;

- Provox Life Home HME/Go HME/Energy HME/Protect HME: pic. 6-8
- Provox Life Night HME: pic. 9-18

### 5.1 Operating instruction

#### 5.1.1 Provox Life Home HME/Go HME/ Energy HME/Protect HME

##### **To attach**

Hold the device and press lightly until the device connects to the attachment device. See fig. 6.

**To speak**

Press the lid to allow speech. See fig. 7.

**To remove**

Hold the device and pull sideways until the device disconnects. See fig. 8.

## 5.1.2 Provox Life Night HME

**To attach**

Hold the device by the side openings between your thumb and middle finger. With your index finger on the top of the device (where the moon symbol is imprinted), press lightly until the device connects to the attachment device. Carefully rotate the device in the attachment device to verify that it has connected correctly. See fig. 9-16.

**To speak**

Occlude the side openings simultaneously to allow speech. See fig. 17.

**To remove**

Hold the device by the side openings between your thumb and middle finger, allowing your index finger to rest on the top of the device (where the moon symbol is imprinted). Pull lightly sideways until the device disconnects. See fig. 9, 10 and 18.

## 6. Device lifetime and disposal

- HMEs are single use devices that are designed to provide pulmonary rehabilitation both day and night, and are to be replaced in line with daily activities.
- The number and type of HMEs needed during a 24 hour period will vary based on the situations encountered and the frequency of involuntary coughing.

- Typically, one HME is used during the night and about two HMEs during the day. The HME may need to be replaced more often if clogged due to mucus secretions or other contamination. HMEs are not reusable and cannot be rinsed in water or any other solution - this washes out the special salt that is essential for HME effectiveness.
- Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

## 7. Compatible products

Provox Life HMEs are compatible with Provox Life attachments.

## 8. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2024 Coloplast A/S. All rights reserved.
Atos Medical AB