 Manufacturer

TRACOE medical GmbH  
Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany  
Tel.: +49 6136 9169-0, fax: +49 6136 9169-200  
info@tracoe.com, www.tracoe.com

**TRACOE**<sup>®</sup>  
*select*

  
aeris.

**CE**  
0 1 2 3

**TRACOE<sup>®</sup> aeris<sup>®</sup> Balloon Dilation Catheter**  
**TRACOE<sup>®</sup> aeris<sup>®</sup> Ballon-Dilatationskatheter**  
**REF 820**

EN	Instructions for Use
DE	Gebrauchsanweisung
FR	Instructions d'utilisation
IT	Istruzioni d'uso
ES	Instrucciones de uso
PT	Instruções de uso
DA	Brugsanvisning
FI	Käyttöohjeet
NO	Bruksanvisning
SV	Bruksanvisning
NL	Gebruiksaanwijzing
CS	Návod k použití
TR	Kullanma talimatı
RU	инструкция по применению
PL	Instrukcja użycia
EL	Οδηγίες χρήσης
KO	사용 설명서
ZH	用于新生儿
AR	دليل استعمال

**DE** / Achtung. Gebrauchsanweisung beachten · **EN** / Warning. Please read Instructions for Use · **FR** / Attention. Respecter le mode d'emploi · **IT** / Attenzione. Consultare le istruzioni per l'uso · **ES** / Atención. Siga las instrucciones de uso · **PT** / Atenção. Ter atenção às instruções de utilização · **DA** / Vigtigt: Følg brugsanvisningen · **FI** / Huomio. Noudata käyttöohjetta · **NO** / OBS. Følg bruksanvisningen · **SV** / OBS! Se bruksanvisningen · **NL** / Let op. Rekening houden met de gebruiksaanwijzing · **CS** / Pozor. Respektujte návod k použití · **TR** / Dikkat. Kullanma talimatına dikkat ediniz · **RU** / Внимание! Соблюдать инструкцию по применению · **PL** / Ostrzeżenie. Zastrzyj do instrukcji używania · **EL** / Προσοχή. Ακολουθήστε τις οδηγίες χρήσης · **KO** / 주의. 사용 설명서를 준수하십시오 · **ZH** / 注意. 注意使用说明书. تنبيه. يجب مراعاة دليل الاستعمال / **AR** ·

**DE** / Artikelnummer · **EN** / Item number · **FR** / Référence · **IT** / Codice prodotto · **ES** / Número de artículo · **PT** / Número do artigo · **DA** / Artikelnummer · **FI** / Tuotenumero · **NO** / Artikkelnummer · **SV** / Katalognummer · **NL** / Artikelnummer · **CS** / Číslo výrobku · **TR** / Ürün numarası · **RU** / Номер артикула · **PL** / Numer katalogowy · **EL** / Αριθμός είδους · **KO** / 품목 번호 · **ZH** / 产品编号 · رقم السلعة / **AR** ·

**DE** / Chargencode · **EN** / Batch code · **FR** / Code de lot · **IT** / Numero di lotto · **ES** / Código de lote · **PT** / Código do lote · **DA** / Batchkode · **FI** / Eräkoodi · **NO** / Chargekode · **SV** / Satsnummer · **NL** / Lotcode · **CS** / Kód šarže · **TR** / Parti kodu · **RU** / Код партии · **PL** / Kod partii · **EL** / Κωδικός παρτίδας · **KO** / 생산 단위 · **ZH** / 批次编码 · رمز الشحن · **AR** ·

**DE** / Hersteller · **EN** / Manufacturer · **FR** / Fabricant · **IT** / Produttore · **ES** / Fabricante · **PT** / Fabricante · **DA** / Producent · **FI** / Valmistaja · **NO** / Produsent · **SV** / Tillverkare · **NL** / Fabrikant · **CS** / Výrobce · **TR** / Üretici · **RU** / Производитель · **PL** / Wytwórca · **EL** / Κατασκευαστής · **KO** / 제조사 · **ZH** / 生产商 · الشركة المصنّعة / **AR** ·

**DE** / Herstellungsdatum · **EN** / Date of manufacture · **FR** / Date de fabrication · **IT** / Data di produzione · **ES** / Fecha de fabricación · **PT** / Data de fabricação · **DA** / Produktionsdato · **FI** / Valmistuspäivä · **NO** / Produksjonsdato · **SV** / Tillverkningsdatum · **NL** / Fabricagedatum · **CS** / Datum výroby · **TR** / Üretim tarihi · **RU** / Дата изготовления · **PL** / Data produkcji · **EL** / Ημερομηνία κατασκευής · **KO** / 제조일 · **ZH** / 生产日期 · تاريخ التصنيع / **AR** ·

**DE** / Verwendbar bis · **EN** / Use by date · **FR** / Utilisable jusqu'au · **IT** / Utilizzare entro il · **ES** / Fecha de caducidad · **PT** / Utilizável até · **DA** / Anvendes inden · **FI** / Käytettävä ennen · **NO** / Kan brukes til · **SV** / Används före – utgångsdatum · **NL** / Bruikbaar tot · **CS** / Použitelné do · **TR** / Son kullanma tarihi · **RU** / Срок годности · **PL** / Użyć do daty · **EL** / Χρήση έως · **KO** / 유통기한 · **ZH** / 保质期至 · قبل للاستخدام حتى / **AR** ·

**DE** / Sterilisiert mit Ethylenoxid · **EN** / Sterilized with Ethylene oxide · **FR** / Stérilisé à l'oxyde d'éthylène · **IT** / Sterilizzato con ossido di etilene · **ES** / Esterilización con óxido de etileno · **PT** / Esterilização com óxido de etileno · **DA** / Sterilisation med ethylenoxid · **FI** / Steriloitu etyleenioksidilla · **NO** / Sterilisering med etylenoksid · **SV** / Sterilisering med etylenoksid · **NL** / Sterilisatie met ethyleenoxide · **CS** / Sterilizováno ethylenoxidem · **TR** / Etilenoksit ile sterilize edilmiştir · **RU** / Стерилизовано оксидом этилена · **PL** / Sterylizowany tlenkiem etylenu · **EL** / Αποστείρωση με αιθυλενοξείδιο · **KO** / 에틸렌으로소독살균되었습니다 · **ZH** / 乙撑氧 (Ethylene oxide) 灭菌 · معقم.مادة أوكسيد الإيثيلين · **AR** ·

**DE** / Nicht erneut sterilisieren · **EN** / Do not resterilize · **FR** / Ne pas restériliser · **IT** / Non risterilizzare · **ES** / No esterilizar de nuevo · **PT** / Não esterilizar de novo · **DA** / Må ikke resteriliseres · **FI** / Ei saa steriloita uudelleen · **NO** / Ikke sterilisert på nytt · **SV** / Får ej återsteriliseras · **NL** / Niet opnieuw steriliseren · **CS** / Znovu nesterilizujte · **TR** / Tekrar sterilize etmeyiniz · **RU** / Не стерилизовать повторно · **PL** / Nie resterylizować · **EL** / Αποστειρωμένο με ακτινοβολία · **KO** / 방사선으로 멸균처리됨 · **ZH** / 经辐射消毒 · معقمة بالإشعاعي / **AR** ·

**DE** / Nicht wiederverwenden · **EN** / Do not reuse · **FR** / Ne pas réutiliser · **IT** / Non riutilizzare · **ES** / No reutilizar · **PT** / Não reutilizar · **DA** / Må ikke genanvendes · **FI** / Ei saa käyttää uudelleen · **NO** / Må ikke benyttes igjen · **SV** / Endast för engångsbruk · **NL** / Niet opnieuw gebruiken · **CS** / Znovu nepoužívejte · **TR** / Tekrar kullanmayınız · **RU** / Не использовать повторно · **PL** / Nie używać powtórnie · **EL** / Μην επαναχρησιμοποιείτε · **KO** / 재사용하지 마십시오 · **ZH** / 不得回收利用 · يُمنع إعادة استخدامها · مرة أخرى / **AR**

**DE** / Trocken aufbewahren · **EN** / Store in a dry place · **FR** / Conserver au sec · **IT** / Conservare in luogo asciutto · **ES** / Conservar en un lugar seco · **PT** / Guardar em local seco · **DA** / Opbevares tørt · **FI** / Säilytettävä kuivassa · **NO** / Oppbevares på et tørt sted · **SV** / Förvaras tørt · **NL** / Droog bewaren · **CS** / Uchovávejte v suchu · **TR** / Kuru şekilde saklayınız · **RU** / Хранить в сухом месте · **PL** / Chronić przed wilgocią · **EL** / Να φυλάσσεται σε ξηρό μέρος · **KO** / 건조한 곳에 보관하십시오 · **ZH** / 需保存于干燥处 · يجب الحفاظ عليها جافة / **AR**

**DE** / Von Sonnenlicht fernhalten · **EN** / Protect from sunlight · **FR** / Protéger de l'ensoleillement · **IT** / Non esporre alla luce del sole · **ES** / Mantener alejado de la luz solar · **PT** / Manter afastado da luz solar · **DA** / Beskyttes mod sollys · **FI** / Säilytettävä auringonvalolta suojattuna · **NO** / Må holdes unna sollys · **SV** / Skyddas från solljus · **NL** / Beschermen tegen zonlicht · **CS** / Chraňte před slunečním světlem · **TR** / Güneş ışığından uzak tutunuz · **RU** / Беречь от солнечных лучей · **PL** / Trzymać z dala od światła słonecznego · **EL** / Να φυλάσσεται μακριά από το φως του ήλιου · **KO** / 직사광선에 노출되지 않도록 하십시오 · **ZH** / 避免阳光直射 · يجب إبقاؤها بعيدة عن أشعة الشمس / **AR**

**DE** / Bei beschädigter Verpackung nicht verwenden · **EN** / Do not use if package is damaged · **FR** / Ne pas utiliser en cas d'emballage endommagé · **IT** / Non utilizzare se la confezione è danneggiata · **ES** / No utilizar si el envase está dañado · **PT** / Não utilizar em caso de embalagem danificada · **DA** / Må ikke anvendes, hvis emballagen er beskadiget · **FI** / Älä käyttää, jos pakkaus on vahingoittunut · **NO** / Må ikke benyttes hvis emballasjen er skadet · **SV** / Används inte om förpackningen är skadad · **NL** / Niet gebruiken bij beschadigde verpakking · **CS** / Nepoužívat při zjištění poškození obalu · **TR** / Ambalajı hasarlıysa kullanmayınız · **RU** / Не использовать, если целостность упаковки нарушена · **PL** / Nie używać, jeżeli opakowanie jest uszkodzone · **EL** / Να μην χρησιμοποιείτε σε περίπτωση που η συσκευασία έχει υποστεί ζημιά · **KO** / 포장이 손상된 경우 사용하지 마십시오 · **ZH** / 如包装破损, 不得使用 · مَنع الاستخدام في حالة تلف العبوة · **AR**

## Instructions for Use

### TRACOE® aeris® Balloon Dilation Catheter

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE.

**Caution:** This device is for operating room use only.

#### Device Description

The TRACOE aeris Balloon Dilation Catheter is comprised of a single lumen catheter with a high pressure balloon near the distal tip. A stylet is provided to facilitate advancement of the balloon dilation catheter to the desired location. The stylet must be removed before inflation of the high pressure balloon. A luer at the proximal end is used for placement of the stylet and injecting sterile water into the balloon. Two radiopaque markers, located on the catheter inside the balloon can be used to confirm balloon placement under Fluoroscopy.

REF	BALLOON SIZE (diameter x length)	MAX INFLATION PRESSURE	INFLATION DEVICE/ SYRINGE VOL
820-05	5 x 30 mm	17 atm	10 ml
820-06	6 x 30 mm	17 atm	10 ml
820-07	7 x 30 mm	17 atm	10 ml
820-08	8 x 30 mm	17 atm	10 ml
820-09	9 x 30 mm	17 atm	10 ml
820-10	10 x 30 mm	17 atm	20 ml
820-12	12 x 40 mm	10 atm	20 ml
820-14	14 x 40 mm	10 atm	20 ml
820-16	16 x 40 mm	10 atm	20 ml
820-18	18 x 40 mm	10 atm	20 ml

1 atm = 1.013125 bar = 101325 Pa

#### Indications for Use

The TRACOE aeris Balloon Dilation Catheter is intended to dilate strictures of the airway.

#### Contraindications

Balloon dilation of the airway is contraindicated in any patient whose degree of respiratory failure would not allow the patient to tolerate the manipulation required to accomplish balloon dilation.

Balloon dilation is contraindicated in the presence of:

- significant active bleeding from the site of the proposed dilation
- and/or presence of a known perforation at the site of proposed dilation
- and/or presence of a known fistula between the tracheobronchial tree and esophagus, mediastinum or pleural space

#### Warnings

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call TRACOE medical GmbH.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the integrity of this device and or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and or cause patient infection or cross-infection, including, but not limited to, the transmission of

infectious disease(s) from one patient to another. Contamination of this device may lead to injury, illness or death of the patient.

EN

### **Precautions**

- Use the catheter prior to the expiration date specified on the package.
- Store in a cool dry location.
- This device should be used only by or under the supervision of a physician thoroughly trained in airway balloon dilation. A thorough understanding of the technical principles, clinical application and risks associated with balloon dilation of the airway is necessary before using this device.
- If resistance is met during the procedure do not advance the catheter without first determining the cause of resistance and taking necessary action.
- After use, dispose of product and packaging in accordance with hospital, administrative and or local government policy.
- Intended for single patient use only, DO NOT REUSE.
- Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.
- Appropriate anesthetic techniques to minimize respiratory effort should be used before this deament, as a result of ongoing vigorous respiration while the balloon is dilated and the airway occluded.
- Careful monitoring of patient's oxygen levels is necessary during balloon dilation. Hypoxia may occur as a result of extended periods of occluding the airway.
- During dilation caution should be used to avoid any interference and possible obstruction of other devices such as tracheostomy tubes and endotracheal tubes.
- Compatibility of this device has not been determined for use within the working channel of a bronchoscope.
- The balloon must be inflated with sterile water. Do not pretest or pre-inflate the balloon.
- The catheter lumen is occluded at the distal end and cannot be used to ventilate.
- Use of a balloon that is too large for the target anatomy and may cause damage to the surrounding anatomy.
- Use of an undersized balloon may result in failure to properly treat the target anatomy. Do not try to force movement of the balloon catheter if balloon is inflated.
- If the balloon catheter migrates during inflation of the balloon do not attempt to advance or retract the balloon without deflating the balloon first.
- Do not advance, retract, or hold the balloon catheter or balloon catheter with stylet against resistance. Tissue damage or tissue trauma may occur as a result.
- Do not advance, retract, or hold the balloon catheter or the balloon catheter with stylet against resistance as damage to the device may occur.
- If at any point during the procedure the balloon does not deflate rupture the balloon with a sharp instrument to allow removal.
- Position balloon catheter properly by using direct visualization. Improperly inflating the balloon in the incorrect location may harm the patient.
- The TRACOE aeris Balloon Dilation Catheter is compatible with the TRACOE high pressure inflation device REF 725.

### **Instructions for Use: Preparation**

Before using inspect the pouch for any breach of the packaging to ensure sterile product and that no damage has occurred to product during shipping.

Visualization of the airway using endoscopy or bronchoscopy, flexible or rigid is recommended in order to determine the location of the stricture and guide placement of the balloon.

Selection of the proper balloon size is critical such that the diameter of the balloon fully inflated does not exceed the expected diameter of a healthy airway. Endoscopic direct visualization or imaging can assist in determining a healthy airway diameter. The cartilage skeleton i.e. tracheal rings should be taken into account in limiting the balloon size. A sharp instrument which can reach the site of dilation should be made available in the event that difficulties arise during deflation of the balloon.

1. Open the sterile package and remove the airway balloon catheter with the stylet in place.
2. If using the stylet confirm that it is locked in place onto the luer at proximal end of the balloon dilation catheter. If not, unlock the stylet, remove and place aside.
3. Prepare the TRACOE high pressure inflation device REF 725.
4. Remove the green protective sheath that covers the balloon.
5. Wipe the outer surface of the balloon and catheter with sterile water soaked gauze pad.

### Placement

1. Locate the stricture using, flexible or rigid, bronchoscopy or endoscopy. Fluoroscopy may be used to confirm balloon placement. Two radiopaque markers are located on the catheter inside the balloon.
2. The TRACOE aeries Balloon may be shaped if desired. The Airway Stylet must be in place and locked for shaping. Kinking may occur without Stylet locked in place.
3. Under endoscopic visualization slowly and gently advance the Airway Balloon to the site of the stricture.
4. Center the balloon portion of the device across the area to be dilated. The end of the balloon should be positioned proximal to the stricture.

**Note:** The balloon on the TRACOE aeries balloon dilation catheter will upon initial inflation present two hubs proximally and distally on the balloon. This is to help secure the balloon in place during inflation. The body of the balloon will increase to create a uniform diameter balloon at recommended inflation pressure. This should be considered an adjunct to the procedures and actions normally taken to prevent slippage in either direction of the intended area to be dilated. Firm control of the balloon catheter should be maintained during the entire procedure.

5. Remove the Airway Stylet, if it has not been removed already, and attach the balloon dilation catheter luer to the connecting tube of the Inflation Device with gauge to monitor balloon pressure.
6. Inflate the balloon with sterile water to desired pressure. Monitor the balloon pressure using the gauge on the inflation device. Visualize endoscopy during inflation of the balloon assessing the diameter, shape and the position of the balloon. Ensure the proximal end of the balloon remains proximal to the stricture throughout inflation.

### Warnings

Do not exceed the Maximum Inflation Pressure stated on the label or in the product table in the Device Description section of this Instructions for Use.

If the balloon moves distally or proximally, out of position, at any time during the procedure do not hold the balloon against any resistance. Deflate the balloon and reposition to desired position. Once the balloon has been repositioned re-inflate.

Holding the balloon against resistance may damage the balloon dilation catheter or anatomy. If the balloon does not deflate rupture the balloon with a sharp instrument to allow removal.

**Note:** If the balloon ruptures at any time during the procedure, identified by a rapid decrease in pressure on the inflation gauge or visually. Carefully, deflate the balloon and remove the balloon dilation catheter.

7. As dilation takes place the pressure readings may fluctuate. Adjust the balloon pressure as necessary to meet the desired target.

### **Withdrawing the Balloon Dilation Catheter**

1. To deflate maintain an endoscopic view of the proximal end of the balloon as vacuum is applied using the inflation/deflation device. Completely deflate the balloon, using vacuum, before pulling back on the balloon dilation catheter. Visually confirm endoscopically the balloon is deflated before removing the balloon dilation catheter.

2. If additional inflations are required, slightly inflate the balloon using fingers gently re-wrap the airway balloon in a clockwise motion, compressing the balloon starting at the distal end of the balloon working back as you apply vacuum. Use the inflation device to apply vacuum.

**Note:** After the balloons first inflation the hubs on the distal and proximal ends of the balloon may not appear upon inflation.

**Note:** Only advance or withdraw the balloon dilation catheter when the balloon is completely deflated and without the green protective sheath covering the balloon. Advancing the balloon partially or fully inflated may cause serious damage to anatomy or the device.

**Note:** Do not rotate the balloon dilation catheter after placing it into the patient.

3. Dispose of the catheter in accordance with accepted hospital guidelines.

### **Possible Complications**

Possible complication that may result from airway balloon dilation are bleeding, perforation, injury to vocal cords, rupture, partial or complete, resulting in pneumomediastinum, pneumothorax, mediastinitis secondary to tracheal dilation, chest pain, laryngospasm, bronchospasm, atelectasis, pulmonary edema, airway obstruction due to edema and hypoxia.

### **Returns**

The return of used products will only be accepted if a decontamination certificate is enclosed. This form is available either directly from TRACOE medical GmbH, or [www.tracoe.com](http://www.tracoe.com).

### **General Terms And Conditions**

The sale, delivery and return of all TRACOE products shall comply with our General Terms and Conditions which are available from TRACOE medical GmbH or [www.tracoe.com](http://www.tracoe.com).

### **Patents**

This device and its use are covered by the following US patent: 7,771,446B2.

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