2022-08
801511
4.0
M

**TRACOE®**

 TRACOE medical GmbH
Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany
Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200
info@tracoe.com, www.tracoe.com

 TRACOE medical GmbH
Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany
Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200
info@tracoe.com, www.tracoe.com

TRACOE®
twist




CE
0 1 2 3

TRACOE® *twist* plus Tracheostomy Tubes











TRACOE® *twist* plus Tracheostomiekänäülen

REF 311, REF 312, REF 313,
REF 314, REF 316, REF 888-316

| | |
|----|----------------------------|
| EN | Instructions for Use |
| DE | Gebrauchsanweisung |
| FR | Instructions d'utilisation |
| IT | Istruzioni d'uso |
| ES | Instrucciones de uso |
| PT | Instruções de uso |
| DA | Brugsanvisning |
| FI | Käyttöohje |
| NO | Bruksanvisning |
| HU | Használati utasítás |
| SV | Bruksanvisning |
| NL | Gebruiksaanwijzing |
| HR | Upute za uporabu |
| SL | Navodila za uporabo |
| CS | Návod k použití |
| RO | Instrucțiuni |
| TR | Kullanma talimatı |
| RU | инструкция по применению |
| PL | Instrukcja używania |
| EL | Οδηγίες χρήσης |
| KO | 사용 설명서 |
| ZH | 用于新生儿 |
| AR | دليل استعمال |

| | REF 311 | REF 312 | REF 313 | REF 314 | REF 316 | REF 888- 316 |
|--|------------|------------|------------|------------|------------|--------------------|
|  | 1 | - | - | - | - | - |
|  | - | 1 | - | - | - | - |
|  | - | - | 1 | - | - | - |
|  | - | - | - | 1 | - | - |
|  | - | - | - | - | 1 | - |
|  | - | - | - | - | - | 1 |
| 9  | - | 1 | - | 1 | - | 1 |
| 9a  | 2 | 1 | 2 | 1 | 2 | 1 |
| 6  | 1 | 1 | 1 | 1 | 1 | 1 |
| 10  | - | 1 | - | 1 | - | 1 |
| 11  | - | - | - | - | 1 | 1 |
| 12  | 1 | 1 | 1 | 1 | 1 | 1 |

DE / Achtung. Gebrauchsanweisung beachten · **EN** / Warning. Please read Instructions for Use · **FR** / Attention. Respecter le mode d'emploi · **IT** / Attenzione. Consultare le istruzioni per l'uso · **ES** / Atención. Siga las instrucciones de uso · **PT** / Atenção. Ter atenção às instruções de utilização · **DA** / Vigtigt: Følg brugsanvisningen · **FI** / Huomio. Noudata käyttöohjetta · **NO** / OBS. Følg bruksanvisningen · **HU** / Figyelem. Olvassa el a használati útmutatót · **SV** / OBS! Se bruksanvisningen · **NL** / Let op. Rekening houden met de gebruiksaanwijzing · **HR** / Upozorenje. Pridržavajte se uputa za uporabu · **SL** / Pozor. Upoštevajte navodila za uporabo · **CS** / Pozor. Respektujte návod k použití · **RO** / Atenție. Respectați instrucțiunile de utilizare · **TR** / Dikkat. Kullanma talimatına dikkat ediniz · **RU** / Внимание! Соблюдать инструкцию по применению · **PL** / Ostrzeżenie. Zazrzyj do instrukcji używania · **EL** / Προσοχή. Ακολουθήστε τις οδηγίες χρήσης · **KO** / 주의. 사용 설명서를 준수하십시오 · **ZH** / 注意。注意使用说明书 · **AR** / تنبيه. يجب مراعاة دليل الاستعمال · **HE** / אזהרה. עקוב אחר הוראות השימוש

DE / Medizinprodukt · **EN** / Medical Device · **FR** / Dispositif médical · **IT** / Dispositivo medico · **ES** / Dispositivo médico · **PT** / Dispositivo médico · **DA** / Medicinsk udstyr · **FI** / Lääkinnällinen laite · **NO** / Medisinsk enhet · **HU** / Orvosi eszköz · **SV** / Medicinteknisk produkt · **NL** / Medisch hulpmiddel · **HR** / Medicinski proizvod · **SL** / Medicinski pripomoček · **CS** / Zdravotnický prostředek · **RO** / Dispozitiv medical · **TR** / Tıbbi cihaz · **RU** / Медицинское изделие · **PL** / Wyrób medyczny · **EL** / Ιατροτεχνολογικό προϊόν · **KO** / 의료 기기 · **ZH** / 医疗器械 · **AR** / מכשור רפואי · **HE** / جهاز طبي

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Instructions for Use for TRACOE® *twist plus* Tracheostomy Tubes

NOTE:

- Please read the instructions for use carefully. They are part of the described product and must be available at all times. For your own and your patients' safety, please observe the following safety information.

- The illustrations to which the text refers can be found on the (foldout) illustrated pages at the beginning of these instructions. The numbers and letters in brackets refer to the respective illustrations and product components of the tracheostomy tube. The symbols and icons used are explained on pages 1 to 5.

1. General Information

CAUTION: The product may only be used by physicians and other **trained** staff responsible for the patient's care! The attending doctor must determine the correct size and type.

Intended Purpose: This medical device is a tracheostomy tube.

Functional Description: The tracheostomy tube ensures that respired air passes through the tracheostoma.

- The inflated cuff separates the upper from the lower respiratory tract, preventing the flow of air from the lungs into the mouth and vice versa. This ensures that the patient can breathe only through the tube. The cuff is inflated by filling it with air using the inflation line. When inflated, the cuff seals the space between the trachea and the outer cannula wall. This allows efficient ventilation with ventilation equipment and prevents subglottic secretions from entering the lower respiratory tract.
- Fenestrated cannulas allow some of the air needed for speaking to reach the upper respiratory tract through the fenestration. Airway resistance in the upper respiratory tract is also reduced.
- See below for a more detailed functional description.

Maximum Useful Life of the Tracheostomy Tube:

- with cuff - 29 days,
- without cuff - 5 weeks,

from first use (see Chapter 5 "General Precautions"). This maximum period of use also includes any periods when the tracheostomy tubes or inner cannulas are not in use.

Product for Single-Patient Use: Multiple use by one and the same patient during the period of use is permitted.

Product Info Card: This product is supplied with a product info card with two detachable labels. The labels show product details. Keep the info card separately to facilitate reordering. The detachable label can be affixed to the patient's medical records, for example.

2. General Description

The product consists of an outer cannula (1) with a cuff (2) or without a cuff and an inner cannula with a 15 mm connector (3), which is already inserted in the outer cannula and is locked with a locking ring (4). Two replacement inner cannulas (9/9a), each with a 15 mm connector, are also supplied.

The cuff is inflated or deflated using the filling valve with pilot balloon (2a), which can be found on the inflation line (2b).

The outer cannula is made of a radiopaque polyurethane and is attached to a neck flange (5) that can swivel on 2 axes (cardanic suspension).

The perforated obturator (6) is already inserted in the inner cannula and makes it easier to insert the tube. The perforation also allows the use of a Seldinger wire ($\varnothing 1.27 \pm 0.04$ mm, e.g. REF 518) for cannulation/recannulation.

3. Indications

The products are suitable for patients where access to the respiratory tract is required via a tracheostoma.

- Tubes without a cuff are particularly suited for patients who breathe spontaneously but need to use a tube.
- Tubes with a cuff are suitable for patients where access to the respiratory tract is required via a tracheostoma with the trachea sealed.
- Fenestrated tubes facilitate speech in patients with an intact larynx.
- Tubes with a suction line (7) are for use in patients where suction of subglottic secretion is indicated.

4. Contraindications

4.1 Absolute contraindications

- Do not use the occlusion cap (10)/ speaking valve in laryngectomised patients (patients without a larynx) - risk of suffocation!

4.2 Relative contraindications

Relative contraindications are cases where the risks and benefits of the procedure must be carefully evaluated:

- Paediatric use
- Unusually deep-lying trachea (e.g. in obese patients) or difficult anatomical conditions

5. General Precautions

- You are strongly advised to keep a ready-to-use spare tube and several spare inner cannulas at the patient's bedside. Store in a clean and dry state.

- Check the tube for integrity and correct function prior to use/insertion. Check, for example, for the absence of obstruction, cuff seal, correct and stable fit of the inner cannula inside the outer cannula, absence of kinks, stable connection between the tube and neck flange, etc. The cuff material must not be brittle. If the product is damaged, replace with a new product.

- Do not use force on the tracheostomy tube, as this may damage or break it. If the connections at the 15 mm connector are tightly fixed, always use a disconnect wedge approved for tracheostomy tubes.

- Keep the 15 mm connector (3) clean and dry.

- When changing the inner cannula, always ensure that the inflation line (2b) of the cuff is not positioned between the inner and outer cannulas as it may get trapped and damaged.

- When repositioning the patient in bed, ensure that the patient does not lie on the pilot balloon (2a), as this would increase the cuff pressure and may damage the trachea.

- During mechanical ventilation and frequent changes of the patient's position or manipulation of the tube, the inner cannula may become separated from the outer cannula.

- The use of tracheostomy tubes can lead to pressure points and skin irritation. To avoid this, we recommend placing a pad underneath the neck flange (e.g. TRACOE purofoam REF 958 or REF 959 or TRACOE softpad REF 969).

- To avoid damage to the material, the cuff should not be in contact with lidocaine-containing aerosols or any ointments.

- All components of the cuff inflation system must lie freely and without kinks during the cuff pressure check. Otherwise it is possible

that the pressure will be displayed incorrectly.

- Thin cuffs are to some extent permeable to water vapour. Therefore, condensate may accumulate in the cuff. Where quantities are small, this is of no significance. However, if larger quantities of condensate are inadvertently aspirated into the inflation line, it is no longer possible to measure and adjust the cuff pressure correctly and the tube must be replaced. Before removing the tube, remove as much air or water from the cuff as possible using a syringe.
- When using together with other medical devices, the respective instructions for use must be observed. If in doubt, contact the manufacturer.
- The use of cuffed tubes on aircraft is not advisable without permanent cuff pressure regulation, as fluctuations in the cuff pressure may occur during the flight.

6. Warnings

- Do not use if the sterile packaging is not intact.
- Use the locking ring (4) to loosen the inner cannula (C). Do not turn the 15 mm connector (3), as this may cause the inner cannula inside the outer cannula to rotate and the air supply to be interrupted.
- Never use fenestrated inner cannulas for ventilation.
- Use the occlusion cap (10) only with fenestrated tubes (outer + inner cannula) and an uninflated cuff.
- Do not use a speaking valve when the patient is asleep.
- When using laser or electrosurgical treatment equipment, ensure that the tracheostomy tube is at a safe distance. There is a risk of fire and toxic gases and the tube may get damaged.
- Select the correct tube to ensure that the fenestration (if applicable) is positioned in the tracheal lumen at a sufficient distance from the stoma canal. Poor positioning exposes ventilated patients to the risk of emphysema. Granulation tissue may also form or airway resistance may increase if speaking valves or occlusion caps are used.
- If a fenestrated tube is used, there is a risk of increased formation of granulation tissue.
- The cuff pressure can change, for example due to altitude changes (e.g. in an aircraft), or if laughing gas is used as an anaesthetic or during the connection of a handheld manometer.
- Excessive cuff pressure poses a risk of permanent damage to the trachea.
- Excessive cuff pressure may cause cuff herniation.
- Insufficient cuff pressure may result in inadequate ventilation and poses a risk of aspiration.

- Ensure that the correct Luer connectors are used (e.g. **2a + 7a**).
- During insertion and removal of the tube, irritation, the need to cough or bleeding may occur.

7. Adverse Reactions

Pressure points, skin irritation, granulation tissue, the need to cough, difficulty swallowing and bleeding.

8. Inserting the Tube

8.1 Preparing the Tube

1. Check the pack content for completeness (**D**).
2. For tubes with a cuff, check the cuff (**2**) for leakage by inflating it. For this purpose, inflate the cuff with a handheld manometer, for example the TRACOE cuff pressure monitor REF 720, to a pressure of 50 cm H₂O (\approx 36.78 mm Hg) and wait for 1 minute to see whether the cuff becomes deflated. If the cuff is leak-tight, remove all air from the cuff with a syringe. Push the cuff upwards in the direction of the neck flange (**5**) to facilitate insertion. If there is any leakage, use a replacement tube and return the defective product as instructed in Chapter 13.

8.2 Preparing the Patient

- Ensure that the patient is well pre-oxygenated before cannulation or recannulation.
- Slightly overextend the patient's neck, if possible, to facilitate insertion.
- In case complications arise during cannulation or recannulation, a tube with a smaller diameter should be at hand for emergency use. In addition, take appropriate safety measures that allow the doctor to carry out short-term ventilation using translaryngeal intubation or a laryngeal mask.

8.3 Inserting the Tube

Carry out the following steps:

1. When inserting the tube (**1**), take care that the obturator (**6**) is not inadvertently pushed back into the tube. This can easily be prevented by holding both the tube and obturator firmly with one hand. To aid insertion, apply a thin layer of water-soluble lubricant to the protruding part of the obturator at the patient end and to the adjacent section of the tracheostomy tube, including the cuff.
2. Once the tube has been inserted into the patient's trachea, remove the obturator immediately.
3. For ventilation, securely connect the inner cannula with the

standardised 15 mm connector (3) to a ventilation machine.

For Cuffed Tubes:

4. Inflate the cuff via the filling valve (2a). Adjust the cuff pressure to the individual ventilation therapy and check at regular intervals. Typically, the pressure should be between 20 cm H₂O (≈ 15 mm Hg) and 30 cm H₂O (≈ 22 mm Hg).

A fall in cuff pressure occurs when a handheld manometer is attached. This effect is more pronounced with smaller tubes.

5. Attach the neck strap (12) to the neck flange (5) to fix the tube in place on the patient's neck.

6. Check the position of the tube as well as its correct function (see Chapter 6 "Warnings"). If required, carefully suction the trachea and tube to ensure that the airways are unobstructed.

7. Now check the cuff pressure to ensure that the cuff has not been damaged during insertion.

8.4 Removing the Tube

Prepare to remove the tube as follows:

- Bend the patient's head slightly backwards and
- Completely deflate the cuff of cuffed tubes (see Chapter 8.5).

8.5 Deflating the Cuff

Before deflating the cuff, ensure that as little secretion as possible can enter the bronchia.

While the cuff is being deflated, secretions should be removed using a suction catheter inserted through the tube. If a tube with subglottic suction has been used, the subglottic space should also be suctioned before deflating the cuff, see Chapter 9.5.

To deflate the cuff, connect a syringe with a standard male Luer connector to the filling valve (2a). Completely remove all air/any possible condensate

If the cuff cannot be deflated, the following measures should be taken:

- Remove the inner cannula. Ensure that the inflation line (2b) is not kinked and repeat the deflation of the cuff.
- If it is still not possible to remove the air/condensate from the cuff, carefully cut through the inflation line between the 15 mm connector (3) and the neck flange (5).

9. Handling

9.1 Changing the Inner Cannulas

If viscous secretion collects in the inner cannula and cannot be suc-

tioned, thus impeding the airflow, replace the inner cannula with a new or cleaned inner cannula.

Loosen the inner cannula by turning (C) the locking ring (4) and remove. Once a new inner cannula has been inserted into the outer cannula, lock in place by turning (C) the locking ring (4) until it clicks into place.

When inserting the inner cannula, ensure that the inflation line (2b) of the cuff is not lying between the inner and outer cannulas, otherwise it may get trapped and damaged.

9.2 Fenestrated Tracheostomy Tubes

Unfenestrated inner cannulas have white locking rings and are used, for example, during ventilation.

To allow the patient to speak, insert a **fenestrated inner cannula** (blue locking ring) **in the fenestrated outer cannula**. Once the cuff has been deflated, a speaking valve can be placed on the 15 mm connector of this inner cannula. The instructions for use of the respective speaking valve must be observed.

9.3 Occlusion Cap of the Fenestrated Tube for Weaning

The fenestrated inner cannula can be closed by attaching the enclosed occlusion cap (10) to the 15 mm connector.

CAUTION: When preparing the tube, ensure that the patient's upper respiratory tract is unobstructed. The upper respiratory tract can be cleared by coughing out or suctioning the secretion. The cuff must be deflated. When closing the tube, the patient's breathing and vital signs must be monitored by a nursing professional. If any signs of breathing difficulties are detected, the occlusion cap must be removed immediately.

9.4 Keeping the Fenestration Open

If the tube is used over a prolonged period, it must be ensured at regular intervals that the fenestration is not blocked by secretions, incrustation or in-growing tissue. If required, replace the tube.

9.5 Application of REF 316, REF 888-316 with Subglottic Suction Line

These tubes are fitted with a soft smooth flat plastic element (7) that is fixed to the outer curve of the outer cannula and ends with two openings immediately above the cuff (2). A suction line (7b) is inte-

grated into this element. The suction line has a female Luer connector (**7a**) at its free end, through which any secretions can be removed using a syringe. Alternatively, a special suction device with a vacuum controller connected to the enclosed connectors (**11**) can be used. After suctioning, ensure that the Luer connector (**7a**) is closed.

CAUTION:

- Ensure during suctioning that no excessive vacuum is used over a prolonged period (maximum of 200 mbar).
- To prevent the subglottic area from drying out, we recommend the use of intermittent suction.
- The suction line may get blocked due to accumulated secretions or adhesions to the tracheal mucosa. If it is intended **to rinse the subglottic suction line** (e.g. with a small amount of air or physiological saline solution), ensure beforehand that the cuff is sufficiently inflated (risk of aspiration).
- The subglottic suction line may cause pressure points in the area of the tracheostoma or may cause the tracheostoma to become irregular. In such cases, the doctor must decide whether this type of tube can still be used.

10. Cleaning, Disinfection and Storage

10.1 Cleaning

The tube (outer and inner cannulas) must be cleaned before reuse. Lukewarm potable water can be used for cleaning. TRACOE medical recommends the use of the cleaning products offered by TRACOE for cleaning. Once the tubes have been cleaned, rinse with potable water and air-dry. When cleaning the outer cannulas, take care not to damage the cuff.

CAUTION:

- Never use aggressive domestic cleaning agents, cleaning agents for dentures or solvents with a high alcohol content, for example, as these may have a lasting adverse effect on the function.
- Do not heat the tubes to temperatures above 60°C, as the safety of the products can then no longer be guaranteed.

10.2 Disinfection

Should disinfection appear necessary TRACOE medical can provide a list of disinfectants on request. The use of unsuitable disinfectants can result in damage to the product.

Once the tubes have been disinfected, rinse with sterile water, distilled water or potable water and air-dry.

10.3 Storage

On delivery from the manufacturer, store the tracheostomy tubes until use in a dry place, as specified in the instructions on the original packaging.

Store cleaned tubes in a clean closed container, in a clean dry place protected from sunlight.

11. Product Modifications

Modifications to the TRACOE products may only be carried out by employees or agents of TRACOE medical GmbH.

12. Disposal

These products may only be disposed of in accordance with the valid national regulations governing waste materials.

13. Returns

Returned used products will only be accepted subject to prior consultation and if a completed decontamination certificate is enclosed. This form is available either directly from TRACOE medical or via the website www.tracoe.com

14. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall be effected exclusively on the basis of the valid General Terms and Conditions (GTC), which are available either from TRACOE medical GmbH or on the website www.tracoe.com.

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