**TRACOE *percutan* Seldinger Guide Wire
(REF 517, REF 518)**

DE	Gebrauchsanweisung
EN	Instructions for Use
FR	Instructions d'utilisation
IT	Istruzioni d'uso
ES	Instrucciones de uso
PT	Instruções de uso
DA	Brugsanvisning
FI	Käyttöohje
NO	Bruksanvisning
HU	Használati utasítás
SV	Bruksanvisning
NL	Gebruiksaanwijzing
HR	Upute za uporabu
SL	Navodila za uporabo
CS	Návod k použití
RO	Instrucțiuni
TR	Kullanma talimatı
RU	инструкция по применению
PL	Instrukcja użytkowania
EL	Οδηγίες χρήσης
KO	사용 설명서
ZH	使用说明
HE	הוראות לשימוש
AR	دليل استعمال

TRACOE medical GmbH

Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany

Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200

info@tracoe.com, www.tracoe.com

DE / Achtung. Gebrauchsanweisung beachten · **EN** / Warning. Please read Instructions for Use · **FR** / Attention. Respecter le mode d'emploi · **IT** / Attenzione. Consultare le istruzioni per l'uso · **ES** / Atención. Siga las instrucciones de uso · **PT** / Atenção. Ter atenção às instruções de utilização · **DA** / Vigtigt: Følg brugsanvisningen · **FI** / Huomio. Noudata käyttöohjetta · **NO** / OBS. Følg bruksanvisningen · **HU** / Figyelem. Olvassa el a használati útmutatót · **SV** / OBS! Se bruksanvisningen · **NL** / Let op. Rekening houden met de gebruiksaanwijzing · **HR** / Upozorenje. Pridržavajte se uputa za uporabu · **SL** / Pozor. Upoštevajte navodila za uporabo · **CS** / Pozor. Respektujte návod k použití · **RO** / Atenție. Respectați instrucțiunile de utilizare · **TR** / Dikkat. Kullanma talimatına dikkat ediniz · **RU** / Внимание! Соблюдать инструкцию по применению · **PL** / Ostrzeżenie. Zazrzyj do instrukcji używania · **EL** / Προσοχή. Ακολουθήστε τις οδηγίες χρήσης · **KO** / 주의. 사용 설명서를 준수하십시오 · **ZH** / 注意. 注意使用说明书 · **AR** / اذرهرا. عكوب اآر الراروة השימוש · **HE** / תבית. יכב מראה ליל האסמל

MD

DE / Medizinprodukt · **EN** / Medical Device · **FR** / Dispositif médical · **IT** / Dispositivo medico · **ES** / Dispositivo médico · **PT** / Dispositivo médico · **DA** / Medicinsk udstyr · **FI** / Lääkinnällinen laite · **NO** / Medisinsk enhet · **HU** / Orvosi eszköz · **SV** / Medicinteknisk produkt · **NL** / Medisch hulpmiddel · **HR** / Medicinski proizvod · **SL** / Medicinski pripomoček · **CS** / Zdravotnický prostředek · **RO** / Dispozitiv medical · **TR** / Tıbbi cihaz · **RU** / Медицинское изделие · **PL** / Wyrób medyczny · **EL** / Ιατροτεχνολογικό προϊόν · **KO** / 의료 기기 · **ZH** / 医疗器械 · **AR** / מכשור רפואי · **HE** / جهاز طبي

Rx ONLY

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Instructions for Use TRACOE *percutan* Seldinger Guide Wire

Note: Please read the instructions for use carefully. They are part of the described product and must be available at all times. For your patients' and your own safety, please observe the following information.

1. Intended Use and Indications for Use

The TRACOE Seldinger Guide Wire with Guiding Catheter (REF 517) facilitates the reinsertion when changing the TRACOE *twist* Tracheostomy Tubes with minimally traumatic inserter.

The TRACOE Seldinger Guide Wire (REF 518) facilitates the reinsertion when changing TRACOE *twist* plus tracheostomy tubes with minimally traumatic insertion system, TRACOE *vario* and *vario* XL tracheostomy tubes with minimally traumatic insertion system, and all TRACOE tracheostomy tubes with perforated obturator.

Patient Population: The product is intended for adults and adolescents.

Intended User: The product should only be used by medical staff trained in tracheostomy care. Lay persons are no intended users.

Indications for Use: The product is intended for patients with a tracheostomy tube, who require their tracheostomy tube changed using the Seldinger technique.

Single Patient Use and Useful Life: This is a single use disposable product. Therefore, the actual use is limited to a maximum of 1 hour.

2. General Description

Contents of REF 517:

- 1 Seldinger Guide Wire with Guiding Catheter in sterile packaging.
 - The Guide Wire is delivered in a sheath and with straight inserter for moving it forward.
 - The Guiding Catheter has a safety stop at its proximal end.

Contents of REF 518:

- 5 Seldinger Guide Wires in individual sterile packaging.
 - Each Guide Wire is delivered in a sheath and with straight inserter for moving it forward.

Supplementary Products:

Supplementary products are listed in section "Supplementary Products".

3. MRI Safety Information

Please see MRI safety information in the respective IFUs of the tracheostomy tube used. The TRACOE Seldinger Guide Wire with Guiding Catheter is MR unsafe. The TRACOE Seldinger Guide Wire is MR unsafe.

4. General Precautions

- When the product is used together with other medical devices, follow their respective instructions for use. Contact the manufacturer if there are any questions, or if assistance is required.
- Safety precautions must be taken in the event of complications during the described procedures, such that a physician can provide immediate ventilation by alternative airways (e.g. trans laryngeal intubation, laryngeal mask).
- It is strongly recommended that a ready-to-use spare device is available.
- The product should be inspected for integrity and function prior to use. If the product is damaged, it should be replaced with a new product.

5. Warnings

- Do not use this product if the sterile packaging has been compromised or damaged, as this may impair the sterility and integrity of the product.
- Refurbishing (including re-sterilization) is not allowed, this may influence the material and function of the product.
- Modification of TRACOE products is not allowed. TRACOE will not be responsible for modified products.

6. Functional Description

Caution:

- The product should be inspected for integrity and function prior to use. If the product is damaged, it should be replaced with a new product.
- It is strongly recommended that a ready-to-use spare device is available.

6.1 Preparation

This is a sterile device, which can be used in an aseptic environment.

The size of the tube and appropriate length subsequently used is determined by a physician.

Prepare the tracheostomy tube, as described in its instructions for use.

The following functions must be checked immediately prior to use. If the device fails the initial inspection, repeat the procedure with a new device. Do not discard the device, and follow instructions provided in chapter "Returns and Complaints".

1. Inspect the sterile packaging to ensure it is secure, undamaged and all components are present.
2. Open the package and inspect the device for damages prior to use.

6.2 Preparing the Patient

Ensure that the patient is optimally pre-oxygenated immediately before insertion or re-insertion of a tracheostomy tube.

6.3 Changing the Tube

The product should be used as follows, unless otherwise specified in the instructions for use of the selected tracheostomy tube.

Caution:

The TRACOE Seldinger Guide Wire should not be inserted deeper into the trachea than the tracheostomy tube as this may harm the trachea or bronchial tissue. Therefore, the guide wire has markings which indicate how far it should be inserted.

6.3a Using Seldinger Guide Wire with Guiding Catheter (REF 517)

1. Prepare the TRACOE *twist* Tracheostomy Tube with minimally traumatic inserter (see chapter "Supplementary Products") by inserting the Guiding Catheter until the safety stop reaches the tip of the inserter.
2. Remove the inner cannula of the tube in situ.
3. Insert the Seldinger Guide Wire into the trachea through the tracheostomy tube in situ by advancing it with the straight inserter. Its proximal marking should be clearly visible outside the tube. Ensure that straight inserter and sheath are complete and discard them.
4. Remove the outer cannula, while securing the Seldinger Guide Wire in situ.
5. Introduce the Guiding Catheter with inserter and tube over the Guide Wire and forward the tube until the distal end of the tube

reaches skin level.

The proximal marking of the Seldinger guide wire should always be visible outside the Guiding Catheter to ensure correct use and positioning.

6. Advance the tracheostomy tube along the Guide Wire into the trachea until the neck plate is in contact with the skin surface. During insertion it is recommended that the tracheostomy tube and inserter/obturator are held together in one hand.

7. While securing the remaining tracheostomy tube at its neck flange, withdraw the minimally traumatic inserter together with the Seldinger Guide Wire and Guiding Catheter. For details please follow the respective instructions for use of the tracheostomy tube.

Warning: After the procedure the inserter must be checked whether it is still undamaged and complete

8. Proceed with tracheostomy tube use, as described in the respective instructions for use.

6.3b Using Seldinger Guide Wire (REF 518)

1. Prepare the TRACOE *twist plus* or *vario/ vario XL* Tracheostomy Tube with minimally traumatic insertion system or perforated obturator (see chapter "Supplementary Products").

2. Remove the inner cannula of the tube in situ where applicable.

3. Insert the Seldinger Guide Wire into the trachea through the tracheostomy tube in situ by advancing it with the straight inserter. Its proximal marking should be clearly visible outside the tube. Ensure that straight inserter and sheath are complete and discard them.

4. Remove the tracheostomy tube, while securing the Seldinger Guide Wire in situ.

When using the minimally traumatic insertion system:

5a. Introduce the Guiding Catheter with inserter and tube over the Guide Wire and forward the tube until the new tube reaches skin level.

The proximal marking of the Seldinger guide wire should always be visible outside the Guiding Catheter to ensure correct use and positioning.

6. Advance the tracheostomy tube along the Guide Wire into the trachea until the neck plate is in contact with the skin surface. During insertion it is recommended that the tracheostomy tube and inserter/obturator are held together in one hand.

7a. While securing the remaining tracheostomy tube at its neck flange, withdraw the minimal traumatic insertion system together with the Seldinger Guide Wire and Guiding Catheter. For details

please follow the respective instructions for use of the tracheostomy tube.

Caution: On no account should the guiding catheter be removed first.

Caution: After the procedure the inserter must be checked whether it is still undamaged and complete.

8. Proceed with tracheostomy tube use, as described in the respective instructions for use.

When using the perforated obturator:

5b. Introduce the perforated obturator and tube over the Guide Wire and forward the tube until the new tube reaches skin level.

6. Advance the tracheostomy tube along the Guide Wire into the trachea until the neck plate is in contact with the skin surface. During insertion it is recommended that the tracheostomy tube and inserter/obturator are held together in one hand.

7b. While securing the remaining tracheostomy tube at its neck flange, withdraw the minimal traumatic insertion system together with the Seldinger Guide Wire and Guiding Catheter. For details please follow the respective instructions for use of the tracheostomy tube.

8. Proceed with tracheostomy tube use, as described in the respective instructions for use.

7. Care and Cleaning

The product is a single-use disposable product. Therefore, no cleaning is required or allowed.

8. Storage

Store the TRACOE products in their original packaging according to the conditions displayed on the packaging.

9. Packaging

The product is provided sterile (with ethylene oxide) which allows insertion under sterile conditions, e.g. for initial insertion in a surgical suite (operating room).

10. Disposal

Used products are to be disposed of in accordance with the valid national regulations, waste management plans, or clinical procedures governing biohazardous waste materials, e.g. the direct disposal in a tear- and moisture-resistant and tight bag or container, which is routed to the local waste disposal system for contaminated medical products.

For further recommendations, contact your hygiene officer in health facilities, or the local waste management for homecare use.

11. Returns and Complaints

Returned products, that have been used, will only be accepted, if TRACOE has agreed to the return and a completed decontamination certificate and complaint report is enclosed. These forms are available either directly from TRACOE medical, or via the website www.tracoe.com.

If the device is involved in a serious incident, as defined in Medical Device Regulation (EU) 2017/745, you may contact either TRACOE medical (complaints@tracoe.com), and the Competent Authority in the country of use.

12. Supplementary Products

12.1 TRACOE Tracheostomy Tubes for REF 517:

TRACOE Tracheostomy Tubes with minimally traumatic inserter

- TRACOE *twist* tracheostomy tube with minimally traumatic inserter (REF 301-P, REF 302-P, REF 306-P and REF 888-306-P; sizes 07-09 available).

12.2 TRACOE Tracheostomy Tubes for REF 518:

TRACOE Tracheostomy Tubes with minimally traumatic insertion system

- TRACOE *twist* plus Tracheostomy Tube with minimally traumatic insertion system (REF 311- P, REF 312-P, REF 316-P, and REF 888-316-P; sizes 07-10 available).
- TRACOE *vario* Tracheostomy Tube with minimally traumatic insertion system (REF 450-P, REF 451-P, REF 460-P, REF 461-P, REF 470- P, and REF 471-P; sizes 07-09 available)

TRACOE Tracheostomy Tubes with perforated obturator

- TRACOE *twist* Tracheostomy Tube with perforated obturator (REF 301-P, REF 302-P, REF 306-P, and REF 888-306-P; sizes 06-09 available)
- TRACOE *twist* Tracheostomy Tube with perforated obturator (REF 301, REF 302, REF 303, REF 304, REF 306, REF 888-306, and REF 309; only sizes 06-10 available)
- TRACOE *twist* plus Tracheostomy Tube with perforated obturator (REF 311-P, REF 312-P, REF 316-P, and REF 888-316-P; sizes 07-10 available)

- TRACOE *twist* plus Tracheostomy Tube with perforated obturator (REF 311, REF 312, REF 313, RREF 314, REF 316, and REF 888-316; sizes 07-10 available)

13. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall be affected exclusively on the basis of the valid General Terms and Conditions (GTC), which are available either from TRACOE medical GmbH or on our website at www.tracoe.com.

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