


# PROVOX HME Life


## Provox Life Home / Go / Protect HME

1.


2.


3.


## Provox Life Night HME


4.


5.






6.

7.

8.

**Disclaimer**

Atos Medical offers no warranty – neither expressed nor implied – to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## **Contents**

EN - ENGLISH .....	7
ES - ESPAÑOL .....	10

## 1. Descriptive Information

### 1.1 Intended Use

Provox Life HMEs are single-use heat and moisture exchangers for patients breathing through a tracheostoma.

### 2. Contraindications

The devices should only be used in accordance with the Instructions for Use. Patients without the physical, cognitive, or mental ability required to attach, remove or operate the devices themselves, should not use the devices independently and should only use them if they are under sufficient supervision of a clinician or a trained caregiver.

The devices should not be used by patients with a low tidal volume, as the added dead space may cause CO<sub>2</sub> (Carbon dioxide) retention, see 4. Technical data.

### 3. Description of the Device

Provox Life HMEs are single-use devices for pulmonary rehabilitation. They come in different levels of humidification, breathing resistance and filtration that makes them suitable for different situations.

The different Provox Life HMEs are:

**Home:** when taking it easy.

**Go:** when you are out and about.

**Protect:** when you need protection from bacteria, virus, dust and pollen.

**Night:** when sleeping.

Provox Life HMEs are compatible with Provox Life attachments.

### 4. Technical Data

	Pressure Drop at 30 l/min*(Pa)	Pressure Drop at 60 l/min*(Pa)	Pressure Drop at 90 l/min*(Pa)	Moisture loss at VT=1000 ml* (mg/l)
Home	60	200	450	19.5
Go	30	100	230	22.5
Protect**	55	180	320	23
Night	65	210	470	18.5

HME internal volume (Dead space)\* max 13 ml

\* According to ISO 9360

\*\* Bacterial Filtration Efficiency >98%, Virus Filtration Efficiency >98%. Test method adapted from ASTM F2101.

## 5. Warnings

- Do not exert unintentional physical pressure on the device. Unintentional or accidental closing or compression of the device may obstruct breathing.
- Do not use non-compatible devices or non-original devices as it may cause personal injury or damage to the device. Do not dismantle or alter the device.
- Do not re-attach a device clogged with mucus, as it may lead to obstructed breathing.
- Do not attach a device that appears slippery by mucus or other fluids, as it will be hard to attach in a correct position.

## 6. Precautions

- Always test the function of the device prior to use. The lid should immediately return to its open position after releasing it. (Not applicable for Provox Life Night HME).
- Pay attention when the device needs to be removed or attached during the night as the night-time situation could make it more difficult to orient the device in a correct position, and it may also be difficult to detect if contaminated.
- The device shall not be used beyond 24 hours, is not reusable and shall not be rinsed in water or in any other solutions. This will substantially reduce the function of the device and can increase the risk for infections due to bacterial colonization.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device. The device will become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.

## 7. How to Use

### 7.1 Operating Instruction

#### 7.1.1 Provox Life Home/Go/Protect HME

##### To attach

Hold the device and press lightly until the device connects to the attachment device. See fig. 1.

##### To speak

Press the lid to allow speech. See fig. 2.

##### To remove

Hold the device and pull sideways until the device disconnects. See fig. 3.

## **7.1.2 Provox Life Night HME**

### **To attach**

Hold the device by the side openings between your thumb and middle finger. With your index finger on the top of the device (where the moon symbol is imprinted), press lightly until the device connects to the attachment device. Carefully rotate the device in the attachment device to verify that it has connected correctly. See fig. 4-11.

### **To speak**

Occlude the side openings simultaneously to allow speech. See fig. 12.

### **To remove**

Hold the device by the side openings between your thumb and middle finger, allowing your index finger to rest on the top of the device (where the moon symbol is imprinted). Pull lightly sideways until the device disconnects. See fig. 4, 5, and 13.

## **8. Device Lifetime and Disposal**

- HMEs are single use devices that are designed to provide pulmonary rehabilitation both day and night, and are to be replaced in line with daily activities.
- The number and type of HMEs needed during a 24 hour period will vary based on the situations encountered and the frequency of involuntary coughing.
- Typically, one HME is used during the night and about two HMEs during the day. The HME may need to be replaced more often if clogged due to mucus secretions or other contamination. HMEs are not reusable and cannot be rinsed in water or any other solution - this washes out the special salt that is essential for HME effectiveness.
- Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

## **9. Reporting**

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Keep away from sunlight and keep dry; Mantener seco y alejado de la luz solar

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Packaging is recyclable; Envase reciclable

Indicates that the product is in compliance with European legislation for medical devices; Indica que el producto cumple la legislación europea sobre productos sanitarios

UK Conformity Assessed (UKCA) marking; Marcado UKCA (UK Conformity Assessed)

UK Responsible Person  
Atos Medical UK Ltd Tottle Road  
Cartwright House  
Nottingham  
Nottinghamshire NG2 1RT  
England United Kingdom

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
Atos and the Atos logo are trademarks of Coloplast A/S.  
© 2024 Coloplast A/S. All rights reserved.

**Atos**  
[atosmedical.com](https://atosmedical.com)

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Atos Medical AB