


PROVOX Adhesives


Provox FlexiDerm


1.1

1.2

1.3

1.4

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

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	Provox FlexiDerm
Provox XtraFlow HME	X
Provox XtraMoist HME	X

Contents

EN - ENGLISH.....5

ES - ESPAÑOL.....7

1. Intended Use

Provox Adhesives are single use devices intended for laryngectomized patients breathing through a tracheostoma. The devices are attached to the skin around the tracheostoma in order to provide attachment of components of the Provox HME System.

2. Description of the Device

Provox Adhesives are adhesives for everyday use, they are designed to ensure an airtight attachment for the Provox HME system components. The adhesives consist of an adapter and an adhesive part with peel-off liners.

Provox FlexiDerm is a very flexible material and has the strongest adhesive properties. It is a sticky, yet soft and flexible adhesive.

3. Warnings

- Do not use Provox FlexiDerm directly after laryngectomy.
- The use of Provox Adhesive need to be assessed both on an individual basis and reviewed regularly through the course of radiotherapy treatment.
- Reuse will affect performance and may cause transfer of micro-organisms leading to infections.
- Only use Provox Adhesive with compatible HMEs and accessories.

4. Precautions

- The adhesive may irritate the skin. Stop using the adhesive if skin irritation develops and consult your clinician.
- When cleaning the skin from residual glue, prevent particles/fluids from entering the tracheostoma.

5. How to Use

5.1 Prepare

Inspect that the adhesive is not damaged or has tears in the adhesive or around the base plate. Always clean the skin with soap and water and let the skin dry before applying the adhesive.

5.2 Apply

For guidance on how to apply, see figures 1.1-1.4.

5.3 Remove

Peel the adhesive off gently, using the finger lift tab. An adhesive remover, may be helpful for removing adhesives. Always clean the skin with soap and water after removal. Dry the area carefully.

The products are intended to be stored in the box until removed for use, as the box protects the products and the label on the box contains important information like batch code, use-by date and a unique barcode.

6. Disposal

Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

7. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Keep away from sunlight and keep dry; Mantener seco y alejado de la luz solar

Instructions for use; Instrucciones de uso

Medical Device; Producto sanitario

Caution, consult instructions for use; Precaución, consultar instrucciones de uso

Store at room temperature. Temporary deviations within the temperature range (max-min) are allowed; Almacenar a temperatura ambiente. Se permiten desviaciones temporales dentro del rango de temperatura (máx.-mín.)


UK Conformity Assessed (UKCA) marking; Marcado UKCA (UK Conformity Assessed)

UK Responsible Person
Atos Medical UK Ltd Tottle Road
Cartwright House
Nottingham
Nottinghamshire NG2 1RT
England United Kingdom

Indicates that the product is in compliance with European legislation for medical devices; Indica que el producto cumple la legislación europea sobre productos sanitarios;

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2024 Coloplast A/S. All rights reserved.

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Atos Medical AB