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Instructions for Use for TRACOE experc Dilation Set for Percutaneous Tracheostomy (REF 520)

Note: Please read the instructions for use carefully. They are part of the described product and must be available at all times. For your own and your patients' safety, please observe the following safety information.

Fold out pages: The illustrations to which the text refers can be found on the (foldout) illustrated pages at the beginning of these instructions. The numbers indicate product components and refer to the respective illustrations of the product. Symbols and icons used with the product are explained in sections "General Description" and "Functional Description".

1. Intended Use and Indications for Use

TRACOE experc Dilation Set for Percutaneous Tracheostomy is intended for performing a percutaneous dilation tracheostomy, applying the Ciaglia technique with a Seldinger guide wire.

Clinical Benefit: With the TRACOE experc Dilation Set an artificial access to the trachea can be created using percutaneous dilation tracheostomy. This allows the placement of a tracheostomy tube for airway management measures, e.g. mechanical ventilation. Moreover, it prevents long-term intubation damages in the laryngeal region. Compared to an open surgical tracheotomy, a percutaneous dilation tracheostomy can be performed bedside, provided that the indications are observed.

Patient Population: The product is intended for adults and adolescents (≥ 18 years).

The product is intended for patients requiring a tracheostomy tube, and who are not indicated for surgical tracheotomy.

Intended User: Percutaneous dilation tracheostomy with the TRACOE experc Dilation Set for Percutaneous Tracheostomy should only be performed by physicians or under the instruction and supervision of physicians who are trained and experienced in the Ciaglia technique. Bronchoscopic monitoring during this procedure is recommended. This intervention should be performed by at least two

people. One person performs the intervention, while the second person administers the anesthetic, monitors the patient and performs the bronchoscopy.

The TRACOE experc Dilation Set for Percutaneous Tracheostomy is intended to be used under sterile conditions.

Indications for Use: This medical device is indicated for cases where access to the respiratory tract is required by means of a percutaneous dilation tracheostomy, with subsequent placement of a tracheostomy tube.

Single Patient Use and Useful Life: This is a single use disposable product. Therefore, the actual use is limited to a maximum of 1 hour.

2. General Description

The TRACOE experc Dilation Set for Percutaneous Tracheostomy is used to create a dilative tracheostoma for the percutaneous insertion of a tracheostomy tube applying the one-step dilation technique according to the Ciaglia method with the Seldinger guide wire.

The dilator is coated with a hydrophilic gel, which reduces friction, as compared to a non-coated dilator.

It is designed to be used together with TRACOE *twist* Tracheostomy Tubes with minimally traumatic inserter, TRACOE *twist* plus Tracheostomy Tubes with minimally traumatic insertion system, or TRACOE *vario* Tracheostomy Tubes with minimally traumatic insertion system (see section "Supplementary Products").

These tracheostomy tubes are available in different diameters and lengths. The models are provided with the cuff deflated. Their minimally traumatic inserter respective minimally traumatic insertion system is remounted. The appropriate diameter and length of the tube is to be determined by the physician.

Fold out pages: Image 1 shows the content of the TRACOE experc Dilation Set for Percutaneous Tracheostomy

1	scalpel
2	syringe
3	puncture needle (14 gauge) with teflon catheter

4	Seldinger guide wire with sheath and straight inserter
5	short dilator (14 Ch/FR)
6	guiding catheter with safety stop
7	TRACOE experc Dilator
8	4 dressings

Supplementary Products:

- Available combinations with TRACOE Tracheostomy Tubes are listed in section "Supplementary Products".
- Products, which can be used in combination with the described one are listed in section "Supplementary Products".

3. MRI Safety Information

Please see MRI safety information in the respective IFU of the TRACOE Tracheostomy Tube used. The TRACOE experc Dilation Set for Percutaneous Tracheostomy is MR unsafe.

4. Contraindications

4.1 Absolute contraindications:

- positioning the tracheostomy tube in emergencies
- application in pediatrics
- existing infections in the area of the tracheostomy
- existing malignancy in the area of the tracheostomy
- unstable fracture of the cervical spine
- uncertainty about the identification of the anatomical orientation points

4.2 Relative Contraindications:

The following list provides relative contraindications for which the user carefully has to consider the risks against the benefits of the procedure:

- an enlarged thyroid gland
- previous surgical procedures in the neck region (e.g. thyroid-ectomy)
- increased tendency to bleed, e.g. in the case of treatment with anticoagulants
- unusually deep-seated trachea, e.g. in the case of obesity

5. General Precautions

- When the product is used together with other medical devices, follow their respective instructions for use. Contact the

manufacturer if there are any questions, or if assistance is required.

- Safety precautions must be taken in the event of complications during the described procedures, such that a physician can provide immediate ventilation by alternative airways (e.g. trans laryngeal intubation, laryngeal mask).
- It is strongly recommended that a ready-to-use spare device is available near the patient.
- The product should be inspected for integrity and function prior to use/insertion. If the product is damaged, it should be replaced with a new product.
- The sterile packaging should be inspected for damage prior to opening. If the packaging is damaged or has been unintentionally opened, the device should not be used.
- It is recommended to carry out the percutaneous tracheotomy under bronchoscopy control, in order to monitor the intra-tracheal position of the puncture needle, the guide wire, the dilator and the tracheostomy tube and to avoid the danger of para-tracheal insertion or injury to the trachea.
- It is recommended to identify tracheal and pre-tracheal conditions (e.g. thyroid gland isthmus, large vessels, cricoid cartilages) with ultrasound, in order to ascertain the optimal spot for puncture and to minimize bleeding.
- Hemorrhages can cause airway obstructions, if the blood is not immediately suctioned.
- Improper storage conditions may result in product or sterile barrier damage.

6. Warnings

- Do not use this product if the sterile packaging has been compromised/damaged, e.g. open edges, holes in packaging etc.
- Refurbishing (including re-sterilisation) is not allowed, this may influence the material and function of the product. The product is single use only.
- Modifications of TRACOE products are not allowed. TRACOE will not be responsible for modified products.

7. Complications and Side Effects

Based on the published literature regarding this type of device, the following complications/side-effects may occur during a percutaneous dilation tracheostomy: Very common ($\geq 10\%$): minor bleeding / Common (1% - 10%): fracture of tracheal cartilages, oxygen desaturation / Uncommon (0.1% - 1%): severe bleeding, hypotension, increase of CPP, cardiac arrhythmia / Very rare ($< 0.01\%$): pneumothorax, pneumomediastinum procedural-related death / Frequency

not known: bronchospasm, injury of the posterior tracheal wall, emphysema, hemodynamic instability, tracheoesophageal fistula, or tracheostomy-related sepsis.

8. Functional Description

Caution:

- During surgery the patient's breathing must be maintained with an endotracheal tube and continually monitored by oximetry.
- Puncture of the trachea must be checked by the aspiration of air.
- Following the puncture, the tip of the puncture needle should be angled at approx. 60° in relation to the anterior wall of the trachea, in order to avoid injuries to the posterior wall of the trachea. This also applies to the subsequent insertion of the guide wire and dilator.
- To avoid injuries to the posterior wall of the trachea the distal tip of the guiding catheter must never extend over the tip of the guide wire. Therefore, the proximal marking of the guide wire must always be at the proximal end of the guiding catheter.
- The safety stop of the guiding catheter should always be placed right in front of the tip of the experc dilator or the inserter/insertion system. The guiding catheter has markings which provide confirmation of the appropriate insertion depth.
- The experc dilator should not be inserted beyond the "MAXIMUM INSERTION" mark (skin level) in order to avoid injury to the trachea at the level of the carina. This is particularly critical with small patients (height approx. 150 cm or smaller; body weight 50 kg or less).
- When inserting the experc dilator or the tracheostomy tube, excessive rotating movements must be avoided as this can damage the trachea.
- Dilation should be performed to match the chosen tube size, so that insertion of the tracheostomy tube is possible without exerting pressure.

8.1 Preparation

- This is a sterile device, intended for use in an aseptic environment.
- The size and appropriate length of the tube subsequently used should be determined by a physician.
- Prepare the tracheostomy tube, as described in its instructions for use.

The following functions must be checked immediately prior to use. If the device fails the initial inspection, repeat the procedure with a new device. Do not discard the device, and follow instructions provided in section "Returns and Complaints".

1. Inspect the sterile packaging to ensure it is secure, undamaged and all components are present.
2. Open the package and inspect the device for damages prior to use.
3. Verify that the guide wire can be advanced without obstruction through the short dilator (14 Ch/FR) and the guiding catheter.
4. Pour sterile water or saline solution into the available cavity of the blister tray.
5. Immerse the narrow end of the experc dilator into the sterile water or saline to activate the lubricant coating, from the tip to the "MAXIMUM INSERTION" mark. The provided dressings can be used to moisten the anterior part of the experc dilator.
6. For preparation of the air aspiration testing, partially fill the syringe with sterile saline solution or sterile water and attach it to the puncture needle with a teflon catheter.

8.2 Preparing the Patient

Ensure that the patient is optimally pre-oxygenated immediately before insertion or re-insertion.

1. The patient should be examined for anatomical anomalies that may adversely affect this procedure.
2. Patient position: supine with a pillow beneath the shoulders and with the head and neck fully extended (See Fig. 2); the head end of the patient's bed should be raised by 30-40°.
3. Disinfect the skin in the neck region and cover it with sterile drapes.
4. Sedate the patient and monitor in accordance with hospital guidelines.
5. Prior to the procedure, the pharynx should be suctioned. If an endotracheal tube is used for ventilation, it should be unblocked and the tube with the cuff withdrawn as far as the entrance of the *larynx* (approx. 1 cm below). Withdrawing the tube prevents the cuff from being pierced during tracheal puncture. The cuff of the endotracheal tube is re-inflated if necessary and ventilation adjusted with respect to volume per minute.

8.3 Procedure

1. Administer a local anesthetic in the region of interest.
2. Typically, the puncture of the trachea is done between tracheal ring 2 and 3, or 3 and 4, in caudal direction. Therefore, push the puncture needle with syringe attached forward in a posterior caudal direction until air bubbles can be aspirated.
3. The puncture needle should now be tilted at an angle of approx. 60° in relation to the anterior wall of the trachea (see Fig. 5). Confirm the position of the cannula tip in the trachea by aspiration of air through the syringe.
4. Make a 1.5 - 2 cm long incision across the placement site (see Fig. 4) to facilitate the dilation of the skin.
5. Withdraw the puncture needle from the trachea, while the Teflon catheter remains in place. Push the catheter forward several millimeters in a posterior and caudal direction (see Fig. 6).
6. Take the guide wire out of its sheath (approx. 8 cm) and advance it until the J-shaped tip is located in the straight inserter.
7. Put the inserter into the Teflon catheter and introduce the guide wire carefully through the Teflon catheter up to its first marking (approx. 10 cm) into the trachea (see Fig. 7). It should be possible to easily push the guide wire forward and move it freely within the catheter. If the guide wire kinks, replace it immediately with a new one. Check, if possible, the orientation of the guide wire through the bronchoscope.
8. Then, the Teflon catheter can be completely removed by pulling it out over the guide wire, leaving the guide wire in its current position within the trachea.
9. The short dilator (14 Ch/FR) is inserted over the guide wire into the trachea. Dilate the puncture site (see Fig. 8). Ensure that the dilator is aligned with the guide wire (to avoid damaging of the tip) and that it does not damage the posterior wall of the trachea. Monitor by bronchoscopy.
10. Remove the short dilator while securing the guide wire in place.
11. Push the guiding catheter with the safety stop in front (in the direction of the arrow) over the guide wire until the safety stop of the guiding catheter reaches skin level (see Fig. 9). The proximal marking of the guide wire should match the proximal end of the guiding catheter.

Warning: It is recommended not to insert the dilator deeper than the "MAXIMUM INSERTION" mark at skin level, to prevent damage of the trachea and/or carina. Avoid rotational movements during insertion of the dilator. In addition, trauma to the tracheal rings should be avoided.

12. Insert the exprec dilator over the guiding catheter and advance until its distal tip reaches the safety stop. In order to expand the puncture channel, push the dilator and the guiding catheter together carefully forward and backward in the trachea until the opening is slightly larger than the external diameter of the selected tracheostomy tube (see Fig. 10). For orientation purposes, there are 38 Ch/FR and 41 Ch/FR ("MAXIMUM INSERTION") markings on the dilator. Ensure that the position of the guiding catheter and guide wire remain constant in relation to the dilator during dilation.
13. Remove the exprec dilator (See Fig. 11).

Caution: Please comply with the instructions for the selected tracheostomy tube used for the following procedure.

- 14a. If using TRACOE *twist* Tracheostomy Tubes the guiding catheter stays in place together with the guide wire.
- 14b. The TRACOE *twist plus* and *vario* Tracheostomy Tubes include a preassembled guiding catheter. Therefore, only the guide wire stays in place.
15. For insertion of the tracheostomy tube follow the respective instructions for use.

8.4 After Care

The tracheostoma should be covered with a sterile wound dressing as required. Changing the tracheostomy tube within the first few days after the initial tracheostomy should only be performed in exceptional cases using the Seldinger technique and intubation equipment.

9. Care and Cleaning

The product is single-use and disposable, therefore, no cleaning is required or allowed.

10. Storage

Store the TRACOE products in their original packaging according to the conditions displayed on the packaging.

11. Packaging

The product is provided sterile (with ethylene oxide) which allows application under sterile conditions.

12. Disposal

Used products are to be disposed of in accordance with national regulations, waste management plans, or clinical procedures gov-

erning biohazardous waste materials, e.g. the direct disposal in a tear- and moisture-resistant and secure bag or container, which is routed to the local waste disposal system for contaminated medical products.

For further recommendations, contact your hygiene officer in health facilities, or the local waste management for homecare use.

13. Returns and Complaints

Returned products, that have been used, will only be accepted, if TRACOE has agreed to the return in advance. A completed decontamination certificate and complaint report is required. These forms are available either directly from TRACOE medical, or via the website www.tracoe.com.

If the device is involved in a reportable incident, as defined in local medical device legislation, please contact TRACOE medical (complaints@tracoe.com), and the appropriate regulatory body in the country of use.

14. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall be affected exclusively on the basis of the valid General Terms and Conditions (GTC), which are available either from TRACOE medical GmbH or on our website at www.tracoe.com.

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TRACOE Tracheostomy Tubes:

TRACOE *twist* Tracheostomy Tubes

Product Description	Size	Order No.	Piece/ Pack- aging
TRACOE <i>twist</i> tracheostomy tube with low-pressure cuff and minimally traumatic inserter (REF 301-P)	07	REF 301-07-P	1
	08	REF 301-08-P	1
	09	REF 301-09-P	1
TRACOE <i>twist</i> tracheostomy tube fenestrated, with low-pressure cuff and minimally traumatic inserter (REF 302-P)	07	REF 302-07-P	1
	08	REF 302-08-P	1
	09	REF 302-09-P	1
TRACOE <i>twist</i> extract tracheostomy tube with low-pressure cuff, subglottic suction line and minimally traumatic inserter (REF 306-P)	07	REF 306-07-P	1
	08	REF 306-08-P	1
	09	REF 306-09-P	1
TRACOE <i>twist</i> extract tracheostomy tube fenestrated, with low-pressure cuff, subglottic suction line and minimally traumatic insertion system (REF 888-306-P)	07	REF 888-306-07-P	1
	08	REF 888-306-08-P	1
	09	REF 888-306-09-P	1

TRACOE *twist* plus Tracheostomy Tubes

Product Description	Size	Order No.	Piece/ Pack- aging
TRACOE <i>twist</i> plus tracheostomy tube with low-pressure cuff, and minimally traumatic insertion system (REF 311-P)	07	REF 311-07-P	1
	08	REF 311-08-P	1
	09	REF 311-09-P	1
	10	REF 311-10-P	1
TRACOE <i>twist</i> plus tracheostomy tube with low-pressure cuff, double fenestrated, and minimally traumatic insertion system (REF 312-P)	07	REF 312-07-P	1
	08	REF 312-08-P	1
	09	REF 312-09-P	1
	10	REF 312-10-P	1
TRACOE <i>twist</i> plus extract tracheostomy tube with low-pressure cuff, with subglottic suction channel and minimally traumatic insertion system (REF 316-P)	07	REF 316-07-P	1
	08	REF 316-08-P	1
	09	REF 316-09-P	1
	10	REF 316-10-P	1
TRACOE <i>twist</i> plus extract tracheostomy tube with low-pressure cuff, double fenestrated, with subglottic suction channel and minimally traumatic insertion system (REF 888-316-P)	07	REF 888-316-07-P	1
	08	REF 888-316-08-P	1
	09	REF 888-316-09-P	1
	10	REF 888-316-10-P	1

TRACOE *vario* Tracheostomy Tubes

Product Description	Size	Order No.	Piece/ Pack- aging
TRACOE <i>vario</i> tracheostomy tube with low-pressure cuff, spiral reinforcement and minimally traumatic insertion system (REF 450-P)	07	REF 450-07-P	1
	08	REF 450-08-P	1
	09	REF 450-09-P	1
TRACOE <i>vario</i> tracheostomy tube with low-pressure cuff and minimally traumatic insertion system (REF 460-P)	07	REF 460-07-P	1
	08	REF 460-08-P	1
	09	REF 460-09-P	1
TRACOE <i>vario</i> extract tracheostomy tube with low-pressure cuff, subglottic suction line and minimally traumatic insertion system (REF 470-P)	08	REF 470-08-P	1
	09	REF 470-09-P	1

TRACOE *vario* XL Tracheostomy Tubes

Product Description	Size	Order No.	Piece/ Packag- ing
TRACOE <i>vario</i> XL tracheostomy tube with low-pressure cuff, spiral reinforcement and minimally traumatic insertion system (REF 451-P)	07	REF 451-07-P	1
	08	REF 451-08-P	1
	09	REF 451-09-P	1
TRACOE <i>vario</i> XL tracheostomy tube with low-pressure cuff and minimally traumatic insertion system (REF 461-P)	07	REF 461-07-P	1
	08	REF 461-08-P	1
	09	REF 461-09-P	1
TRACOE <i>vario</i> XL extract tracheostomy tube with low-pressure cuff, subglottic suction line and minimally traumatic insertion system (REF 471-P)	08	REF 471-08-P	1
	09	REF 471-09-P	1

TRACOE experc Sets inclusive TRACOE Tracheostomy Tubes:

inclusive TRACOE *twist* Tracheostomy Tubes

Product Description	Size	Order No.	Piece/ Packaging
TRACOE experc Set <i>twist</i> inclusive Dilation Set and <i>twist</i> tracheostomy tube REF 301-P	07	REF 320-07	1
	08	REF 320-08	1
	09	REF 320-09	1
TRACOE experc Set <i>twist</i> inclusive Dilation Set and <i>twist</i> tracheostomy tube REF 302-P	07	REF 321-07	1
	08	REF 321-08	1
	09	REF 321-09	1
TRACOE experc Set <i>twist</i> extract inclusive Dilation Set and <i>twist</i> extract tracheostomy tube REF 306-P	07	REF 322-07	1
	08	REF 322-08	1
	09	REF 322-09	1
TRACOE experc Set <i>twist</i> extract inclusive Dilation Set and <i>twist</i> extract tracheostomy tube REF 888-306-P	07	REF 888-322-07	1
	08	REF 888-322-08	1
	09	REF 888-322-09	1

inclusive TRACOE *twist* plus Tracheostomy Tubes

Product Description	Size	Order No.	Piece/ Packaging
TRACOE experc Set <i>twist</i> plus inclusive Dilation Set and <i>twist</i> plus tracheostomy tube REF 311-P	07	REF 330-07	1
	08	REF 330-08	1
	09	REF 330-09	1
	10	REF 330-10	1
TRACOE experc Set <i>twist</i> plus inclusive Dilation Set and <i>twist</i> plus tracheostomy tube REF 312-P	07	REF 331-07	1
	08	REF 331-08	1
	09	REF 331-09	1
	10	REF 331-10	1
TRACOE experc Set <i>twist</i> plus extract inclusive Dilation Set and <i>twist</i> plus extract tracheostomy tube REF 316-P	07	REF 332-07	1
	08	REF 332-08	1
	09	REF 332-09	1
TRACOE experc Set <i>twist</i> plus extract inclusive Dilation Set and <i>twist</i> plus extract tracheostomy tube REF 888-316-P	07	REF 888-332-07	1
	08	REF 888-332-08	1
	09	REF 888-332-09	1
	10	REF 888-332-10	1

inclusive TRACOE *vario* Tracheostomy Tubes

Product Description	Size	Order No.	Piece/ Pack- aging
TRACOE <i>experc</i> Set <i>vario</i> inclusive Dilation Set and <i>vario</i> tracheos- tomy tube REF 450-P	07	REF 420-07	1
	08	REF 420-08	1
	09	REF 420-09	1
TRACOE <i>experc</i> Set <i>vario</i> inclusive Dilation Set and <i>vario</i> tracheos- tomy tube REF 460-P	07	REF 423-07	1
	08	REF 423-08	1
	09	REF 423-09	1
TRACOE <i>experc</i> Set <i>vario</i> extract inclusive Dilation Set and <i>vario</i> tracheostomy tube REF 470-P	08	REF 422-08	1
	09	REF 422-09	1

 inclusive TRACOE *vario* XL Tracheostomy Tubes

Product Descrip- tion	Size	Order No.	Piece/ Pack- aging
TRACOE <i>experc</i> Set <i>vario</i> XL inclusive Dila- tion Set and <i>vario</i> XL tracheostomy tube REF 451-P	07	REF 421-07	1
	08	REF 421-08	1
	09	REF 421-09	1
TRACOE <i>experc</i> Set <i>vario</i> XL inclusive Dila- tion Set and <i>vario</i> XL tracheostomy tube REF 461-P	07	REF 424-07	1
	08	REF 424-08	1
	09	REF 424-09	1
TRACOE <i>experc</i> Set <i>vario</i> XL extract inclusive Dilation Set and <i>vario</i> XL tracheo- stomy tube REF 471-P	08	REF 425-08	1
	09	REF 425-09	1