


PROVOX Flush

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


1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30					


x 12

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

+

⑫


70% Ethanol ⌚ ·10 min
70% 乙醇 ⌚ 10 分钟

or / 或

70% Isopropylalcohol ⌚ 10 min
70% 异丙醇 ⌚ 10 分钟

or / 或

3% Hydrogenperoxide ⌚ 60 min
3% 过氧化氢 ⌚ 60 分钟

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Due to local Chinese requirements, the text in the Intended use paragraph and the Overall description and product composition paragraph are not translated verbatim.

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Intended Use

Provox Flush is intended to be used to flush drinking water or air through the inner lumen of a Provox voice prosthesis for cleaning purposes. The Provox Flush is intended for both home and clinical use by patient or clinician.

Overall Description and Product Composition (结构及组成)

Provox Voice Prosthesis is composed of Silicone, Polyvinylidene Fluoride (PVDF) and Silicone adhesive. The insertion system in the Provox Vega Puncture Set is composed of MethylMethacrylate Acrylonitrile Butadiene Styrene (MABS), stainless steel, Thermoplastic Styrene-Ethylene/Butylene-Styrene (TPS-SEBS), polypropylene (PP), Polyvinylidene Fluoride (PVDF). The insertion system in the Provox Vega and Provox Vega XtraSeal is composed of Polypropylene (PP), colorant, and Silicone oil.

Provox Flush in the Provox Accessories is composed of Silicone, Polypropylene (PP) and blue masterbatch; Provox Vega Plug in the Provox Accessories is made of Silicone. Provox Voice Prosthesis and its set are sterilized by ethylene oxide. The product is for single use.

The Provox accessories are non-sterile products. The shelf life of package sterilization is 5 years.

Contraindications

None.

Warnings

- The flush is only recommended for patients with sufficient manual dexterity, acceptable vision and satisfactory cognitive ability.

Description of the Device

- Provox Flush is intended to be used to flush drinking water or air through the inner lumen of a Provox voice prosthesis for cleaning purposes.
- The Provox Flush is intended for both home and clinical use by patient or clinician.

Precaution

- The device is intended for single patient use only. Reuse in another patient may cause cross contamination.

Instructions for Use

For instructions how to use the device, see fig. 1-8.

Cleaning and Disinfection

Clean the Provox Flush after each use according to fig. 9-11. Do not use any water other than drinking water to clean and rinse the device.

Disinfect the Provox Flush at least once a day (fig. 12) with one of the following methods:

- Ethanol 70% for 10 minutes
- Isopropylalcohol 70% for 10 minutes
- Hydrogenperoxide 3% for 60 minutes

Caution: Do not clean or disinfect by any other method as this might cause product damage and injury.

If the Provox Flush looks dirty or has air dried in an area with a risk of contamination the device should be both cleaned and disinfected before use. A risk of contamination could exist if the device has been dropped on the floor or if it has been in contact with a pet, someone with a respiratory infection, or any other gross contamination.

During hospitalization it is important to both clean and disinfect the Provox Flush after use but also **before** use since there is an increased risk of infection and contamination. In a hospital it is better to use sterile water for cleaning and rinsing, rather than drinking water.

Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

Product Model

Model
8109-18

Storage Conditions

Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.

Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

Legal Agent and After-sales Service Information

Coloplast (China) Medical Devices Ltd.
Address: Unit1301-1306, 13th Floor, Building 1, No.5
Lido Huayuan Road, Chaoyang District,
Beijing

Agent contact information:

Tel: 010-5920 1888

Fax: 010-5920 1898

Coloplast Customer Service Hotline

Hotline: 400 700 7668

Website: www.coloplast.com.cn

Production Date and Validity

Please see the label for the production date; the shelf life of product is 5 years.

Registration Certificate Number of Medical Devices

Certificate number of medical devices: 国械注进
20223130609

Technical Requirements for Medical Devices

Technical requirement number for medical devices: 国械
注进20223130609

Edition of Instruction for Use

Article number: 11703, Edition date: 2024-06-04

Registration and Manufacturer Company

Atos Medical AB

Registration & Manufacturer location: Kraftgatan 8,
SE-242 35 Hörby, Sweden

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Product reference number; 产品编号

Non-sterile; 未灭菌

Keep away from sunlight and keep dry; 怕雨, 怕晒

Storage temperature limit; 储存温度限制

Store at room temperature. Temporary deviations within the temperature range (max-min) are allowed;
室温下存放。允许温度范围内(最高-最低)的暂时偏差。

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
UK Conformity Assessed (UKCA) marking;
英国合格评定 (UKCA) 标志;

UK Responsible Person
Atos Medical UK Ltd Tottle Road
Cartwright House
Nottingham
Nottinghamshire NG2 1RT
England United Kingdom

Indicates that the product is in compliance with European legislation for medical devices; 法国 Triman 标志和分类信息

Triman symbol and Infotri for France;
法国 Triman 标志和分类信息;

Recycling guidelines; 回收指南;

XXXXX, NN
YYYY-MM-DD

XXXXX, NN = Reference number, Version number;
参考编号, 版本号

YYYY-MM-DD = Date of issue; 发布日期

THIS PAGE IS DELIBERATELY LEFT EMPTY

THIS PAGE IS DELIBERATELY LEFT EMPTY

THIS PAGE IS DELIBERATELY LEFT EMPTY

11703, 03
2024-06-24

Atos
atosmedical.com

Atos Medical AB
Kraftgatan 8, SE-242 35 Hörby, Sweden
Tel: +46 (0)415 198 00 • info@atosmedical.com

Atos and the Atos logo are trademarks of Coloplast A/S.
© 2024 Coloplast A/S. All rights reserved.