
# PROVOX Adhesives

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## 1. StabiliBase 深凸稳固型


Figure 1.1


Figure 1.2


Figure 1.3


Figure 1.4


Figure 1.5


Figure 1.6


Figure 1.7


Figure 1.8

## 2. FlexiDerm 灵活型, XtraBase 微凸加固型, OptiDerm 柔肤型


Figure 2.1


Figure 2.2


Figure 2.3


Figure 2.4

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Due to local Chinese requirements, the text in the Intended use paragraph and the Overall description and product composition paragraph are not translated verbatim.

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ZH - ZHCN..... 11

## 1. Intended Use

Provox Adhesives are single use devices intended for laryngectomized patients breathing through a tracheostoma. The devices are attached to the skin around the tracheostoma in order to provide attachment of components of the Provox HME System.

## 2. Overall Description and Product Composition

Provox Adhesives consist of a baseplate with an adaptor, adhesive tape and peeloff liner. The devices are made of polymer materials and are non-sterile. The devices can be used together with Provox HMEs, Provox Shower Aid or Provox Baseplate Adaptor. When used, the adhesive is attached on intact skin around the tracheostoma. After application a Provox HME, Provox Shower Aid or a Provox Baseplate Adaptor can be attached to the adhesive adaptor, see instruction for use section.

## 3. Description of the Device

Provox Adhesives are adhesives for everyday use, they are designed to ensure an airtight attachment for the Provox HME system components. The adhesives consist of an adapter and an adhesive part with peel-off liners.

**Provox FlexiDerm** is a soft, flexible and acrylic based adhesive with strong adhesive properties.

**Provox XtraBase** is an acrylic based adhesive with strong adhesive properties. It has a concave base which can be especially suitable for somewhat deep tracheostomas.

**Provox StabiliBase** is an acrylic based adhesive with strong adhesive properties. It has a firm base with vertical stabilizing bars that provides support to the tracheostoma during speech. The design of the StabiliBase adapter can be especially suitable for deep tracheostomas. The peel-off liner is divided into three sections to facilitate application to the skin.

**Provox OptiDerm** is made of a hydrocolloid material that is suitable for sensitive skin. For example, dependent on clinical suitability and continued tolerance, it can be used postoperatively and during and after radiotherapy. The hydrocolloid material is commonly used in wound care and one of its main characteristics is that it absorbs wound fluids. When the material has absorbed fluids it becomes white and depending on the amount of secretions it can become jelly-like and it will no longer stick to the skin. Please note that the material may also absorb mucus coming from the stoma, and therefore it remains very important to always wipe and clear the mucus away from your adhesive and stoma area after you cough.

## 4. Contraindications

None.

## 5. Warnings

- Do not use Provox FlexiDerm, Provox StabiliBase or Provox XtraBase directly after laryngectomy.
- The use of Provox Adhesives need to be assessed both on an individual basis and reviewed regularly through the course of radiotherapy treatment.
- Reuse will affect performance and may cause transfer of micro-organisms leading to infections.
- Only use Provox Adhesives with compatible HMEs and accessories.

## 6. Precautions

- The adhesive may irritate the skin. Stop using the adhesive if skin irritation develops and consult your clinician.
- When cleaning the skin from residual glue, prevent particles/fluids from entering the tracheostoma.
- When using OptiDerm in the postoperative period or on sensitive skin, it should be removed very slowly and carefully.

## 7. How to Use

### Prepare

Inspect that the adhesive is not damaged or has tears in the adhesive or around the base plate. Always clean the skin with a Provox Cleaning Towel and/or soap and water and let the skin dry before applying the adhesive.

### Apply

If needed, you can apply skin protection products, like Provox Skin Barrier or adhesive supporting products such as Provox Silicone Glue (read the Instructions for Use accompanying the product). For guidance on how to apply, see figures 1.1 - 1.8 and 2.1 - 2.4.

### Remove

Peel the adhesive off gently, using the finger lift tab. An adhesive remover, e.g. Provox Adhesive Remover (read the Instructions for Use accompanying the product), may be helpful for removing adhesives. Always clean the skin with a Provox Cleaning Towel and/or soap and water after removal. Dry the area carefully.

## 8. Product model

Model	Specification
7253-18	Provox Adhesive FlexiDerm Round (20 pcs)
7254-18	Provox Adhesive FlexiDerm Oval (20 pcs)
7331-18	Provox Adhesive FlexiDerm Plus (20 pcs)
7255-18	Provox Adhesive OptiDerm Round (20 pcs)
7256-18	Provox Adhesive OptiDerm Oval (20 pcs)
7332-18	Provox Adhesive OptiDerm Plus (20 pcs)
7289-18	Provox StabiliBase (15 pcs)
7265-18	Provox XtraBase Adhesive (20 pcs)

## 9. Storage conditions

Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.

## 10. Manufacturer and filing company

Filing & Manufacturer company: Atos Medical AB

Filing & Manufacturer location: Kraftgatan 8, SE-242 35 Hörby, Sweden

Filing & Manufacturer contact: +46 (0415 198 00 • info@atosmedical.com

Production location: Kraftgatan 8, SE-242 35 Hörby, Sweden

Country of origin: Sweden

## 11. Edition of Instruction for Use

Article number: 11617, Edition date: 2024-09-04.

## 12. Disposal

Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

### **13. Reporting**

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

### **14. Legal Agent and After-sales service information**

Coloplast (China) Medical Devices Ltd.  
Address: Unit1301-1306, 13th Floor, Building 1, No.5  
Lido Huayuan Road, Chaoyang District,  
Beijing  
Agent contact information:  
Tel: 010-5920 1888  
Fax: 010-5920 1898  
Coloplast Customer Service Hotline  
Hotline: 400 700 7668  
Website: [www.coloplast.com.cn](http://www.coloplast.com.cn)

### **15. Production date and validity**

Please see the label for the production date; the shelf life of product is 3 years.

### **16. Filing certificate number of medical devices**

Filing certificate number of medical devices: 国械备 20230011

### **17. Technical requirements for medical devices**

Technical requirement number for medical devices: 国械备 20230011

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Keep away from sunlight and keep dry; 怕晒,怕雨

Storage temperature limit; 储存温度限制

Store at room temperature. Temporary deviations within the temperature range (max-min) are allowed;  
室温下存放。允许温度范围内(最高-最低)的暂时偏差。

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Instructions for use; 使用说明书

Medical Device; 医疗器械

Indicates that the product is in compliance with European legislation for medical devices; 表示本产品符合欧洲医疗器械法规

Recycling guidelines; 回收指南

**XXXXX, NN**  
**YYYY-MM-DD**

XXXXX, NN = Reference number, Version number; 参考编号, 版本号


YYYY-MM-DD = Date of issue; 發行日期; 发布日期

**UK\***  
**CA**

UK Conformity Assessed (UKCA) marking; 英國合格評定 (UKCA) 標誌; 英国合格评定 (UKCA) 标志

UK Responsible Person  
Atos Medical UK Ltd Tottle Road  
Cartwright House  
Nottingham  
Nottinghamshire NG2 1RT  
England United Kingdom

\* Not applicable for China

Triman symbol and Infotri for France; 法国 Triman 标志和分类信息;

\* Not applicable for China

11617, 01  
2024-09-04

**Atos**  
atosmedical.com

Atos Medical AB  
Kraftgatan 8, SE-242 35 Hörby, Sweden  
Tel: +46 (0)415 198 00 • [info@atosmedical.com](mailto:info@atosmedical.com)

Atos and the Atos logo are trademarks of Coloplast A/S.  
© 2024 Coloplast A/S. All rights reserved.