
PROVOX

Life

Oxygen

Provox Life Oxygen


Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6


Figure 7


Figure 8

Prescription Information

CAUTION: United States Federal law restricts this device to sale, distribution, and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty—neither expressed nor implied—to the purchaser here under as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Compatible Products; Kompatible Produkte; Compatibele producten; Produits compatibles; Prodotti compatibili; Productos compatibles; Produtos compatíveis; Kompatibla produkter; Kompatible produkter; Kompatible produkter; Yhteensopivat tuotteet; Samhæfar vörur; Ühilduvad tooted; Kompatibilní produkty; Kompatibilis termékek; Zgodne produkty; Kompatibilni proizvodi; Συμβατά προϊόντα; 互換性のある製品; المنتجات المتوافقة

Provox Life HMEs (Home, Go, Energy, Protect, Night)
Provox Life Adhesives (Standard, Sensitive, Night, Stability)
Provox Life LaryTubes
Provox Life LaryButtons

Contents

EN - ENGLISH	7
DE - DEUTSCH	10
NL - NEDERLANDS	14
FR - FRANÇAIS	18
IT - ITALIANO	22
ES - ESPAÑOL	26
PT - PORTUGUÊS	30
SV - SVENSKA	34
DA - DANSK	37
NO - NORSK	40
FI - SUOMI	43
IS - ÍSLENSKA	47
ET - EESTI	50
CS - ČESKY	53
HU - MAGYAR	56
PL - POLSKI	59
HR - HRVATSKI	63
EL - ΕΛΛΗΝΙΚΑ	66
AR - العربية	72

1. Intended Use

Provox Life Oxygen is a single-use accessory to Provox Life HMEs. It allows additional oxygen (O₂) supply for patients breathing through a tracheostoma using Provox Life HMEs.

Intended User Group

Provox Life Oxygen is a device for home and hospital use for adults (18+ years) breathing through a tracheostoma.

2. Contraindications

The device shall not be used by patients with reduced mental or physical cognitive ability. Patients who are unable to attach or remove the device themselves, or without sufficient knowledge how to use the device, or the cognitive ability to understand the risks connected to the use, should not use the device.

The device shall not be used by patients who have a low tidal volume, as the extra dead space may cause carbon dioxide retention.

3. Description of the Device

Provox Life Oxygen is a single-use accessory to Provox Life HMEs for adult patients requiring supplemental oxygen. The device is used between the HME and the Provox Life attachments (see Provox Life Oxygen Overview). Supplemental oxygen (non-humidified) can then be supplied via the oxygen connector of Provox Life Oxygen (see Position 1).

Oxygen tubing with an inner diameter of 3.2 mm (1/8 in.) is recommended. Other diameters of oxygen tubing can be used but must be assessed on a case-by-case basis by the responsible healthcare professional.


Flow rates of up to 6 l/min can be used. However, the HME effect will be reduced as the oxygen flow increases.

Oxygen preserved during inhalation (oxygen reaching the patient) is typically more than ≥78%*.

Note: Use only with dry oxygen.

**Calculated in the worst-case scenario—Provox Life Energy attached to the device and test conditions: tidal volume 250 ml, breathing frequency: 20/min, oxygen flow rate: 6 l/min.*

Provox Life Oxygen Overview

- 1) Oxygen connector
- 2) Compatible with the Provox Life HME range
- 3) Compatible with Provox Life attachments; Adhesives, LaryTubes, LaryButtons

4. Technical Data

	Pressure Drop at 30 l/min of HME + device mounted* (Pa)	Moisture loss at VT=500 ml with supply of O ₂ 3 l/min* (mg/l)	Moisture loss at VT=500 ml with supply of O ₂ 6 l/min* (mg/l)
Home	81	18.3	19.7
Go	32	20.2	21.1
Energy	16	20.9	21.5
Protect	61	21.4	22.6
Night	78	17.1	18.5

Provox Life HME internal volume (dead space) max. 13 ml*

Internal volume (dead space) of Provox Life Oxygen/device with HME attached max. 16 ml

** According to ISO 9360*

5. Warnings

- Do not reuse the device on another patient. Reuse in another patient may cause cross-contamination which may lead to serious infection.
- A high secretion load may indicate more frequent replacement of the device to avoid higher breathing resistance and blocking.

6. Precautions

- Before use, check the device for deformity. Do not use if the device shows any signs of defect. Small parts of the device may enter the airway which may lead to coughing, airway irritation, and possibly pulmonary infection.
- Single use only. The device shall not be used beyond 24 h, is not reusable, cannot be sterilized, and shall not be rinsed in water or in any other solutions. This can increase the risk for serious infection.

Side-Effects

Possible side-effects relating to the use of Provox Life Oxygen may include breathing resistance, coughing, irritation/discomfort of the trachea and/or of the peristomal area.

7. How to Use

Attach

Provox Life HME

1. Hold the device and attach it to the HME (see Figure 1) according to the HMEs own Instructions for Use (IFU).

Oxygen tubing

1. Attach the oxygen tubing securely to the oxygen connector (see Figure 2). Make sure it is pressed all the way.

Note: Attaching the oxygen tubing while the device is in place is more difficult and it can cause unintended movement of the device it is connected to, which may cause coughing, irritation, and discomfort.

To Provox Life attachment

1. Hold the device and press lightly until the device connects to the Provox Life attachment (see Figure 3).
2. Verify that the device is properly secured to the Provox Life attachment (see Figure 4).

Speak

With oxygen tubing connected

1. Occlude the HME according to its own IFU (see Figure 5).

With no oxygen tubing connected

1. Occlude the oxygen connector with a finger and then occlude the HME according to its own IFU (see Figure 6).

Remove

Remove from Provox Life attachment

1. Hold the device and pull sideways until the device disconnects (see Figure 7).

Oxygen tubing

1. Remove the oxygen tubing from the oxygen connector.

Note: Detaching the oxygen tubing while the device is in place is more difficult and it can cause unintended movement of the device it is connected to, which may cause coughing, irritation, and discomfort.

Provox Life HME

1. Hold the device and remove the HME according to its own IFU (see Figure 8).

If further explanation is needed, contact your clinician.

8. Device Lifetime and Disposal

Provox Life Oxygen is for single use and must be replaced at least every 24 h, or more often as needed, e.g. if it is clogged by mucus, becomes dirty, contaminated, or shows any signs of damage. Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

9. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2025 Coloplast A/S. All rights reserved.
Atos Medical AB