
# PROVOX<sup>®</sup> Vega<sup>™</sup>

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



Figure. 19 图19

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Due to local Chinese requirements, the text in the Intended use paragraph and the Overall description and product composition paragraph are not translated verbatim.

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**Provox® Vega™****1. Descriptive information****1.1 Intended use**

Provox Vega Voice Prosthesis is a sterile single use indwelling voice prosthesis intended for voice rehabilitation after surgical removal of the larynx (laryngectomy). Cleaning of the voice prosthesis is performed by the patient while it remains in situ.

The Provox Insertion System is a sterile single use device intended for anterograde replacement of the Provox Vega Voice Prosthesis. This replacement procedure is carried out by a medical doctor or a trained medical professional in accordance with local or national guidelines.

*The Provox Insertion System is not intended to be used for insertion of a voice prosthesis in a freshly made puncture.*

**1.2 Overall description and product composition (结构及组成)**

Provox Voice Prosthesis is composed of Silicone, Polyvinylidene Fluoride (PVDF) and Silicone adhesive.

The insertion system in the Provox Vega Puncture Set is composed of Methyl Methacrylate Acrylonitrile Butadiene Styrene (MABS), Stainless steel, Polyamide (PA66) and epoxy glue, Thermoplastic Styrene-Ethylene/Butylene-Styrene (TPS-SEBS), polypropylene (PP), Polyvinylidene Fluoride (PVDF).

The insertion system in the Provox Vega and Provox Vega XtraSeal is composed of Polypropylene (PP), colorant, and Silicone oil.

Provox Flush in the Provox Accessories is composed of Silicone, Polypropylene (PP) and blue masterbatch; Provox Vega Plug in the Provox Accessories is made of Silicone. Provox Voice Prosthesis and its set are sterilized by ethylene oxide. The product is for single use.

The Provox accessories are non-sterile products.

The shelf life of package sterilization is 5 years.

**1.3 Description of the device General**

Provox Vega is a one-way valve (prosthesis) that keeps a TE-puncture open for speech, while reducing the risk of fluids and food entering the trachea. The Provox Vega voice prosthesis is not a permanent implant, and needs periodic replacement. The prosthesis (Fig. 1) is available in different diameters and several lengths. The device is made of medical grade silicone rubber and fluoroplastic.

**Provox Vega package**

The Provox Vega package contains the following items:

- 1 Provox Vega voice prosthesis pre-assembled in a single-use Insertion System, sterile (Fig. 2)
- 1 Provox Brush of a size corresponding to the voice prosthesis, non-sterile
- 1 Provox Vega Clinician's manual
- 1 Provox Vega Patient's manual
- 1 Provox Brush Instructions for Use

**1.4 CONTRAINDICATIONS**

There are no known contraindications for use or replacement of the Provox Vega voice prosthesis among patients already using prosthetic voice rehabilitation.

**1.5 WARNINGS**

- **Dislodgement or extrusion of the Provox Vega voice prosthesis** from the TE puncture and subsequent ingestion, aspiration or tissue damage may occur. A foreign body in the airway may cause severe complications such as acute respiratory distress and/or respiratory arrest.
- **Select the proper prosthesis size.** A tight fit may cause tissue necrosis and extrusion.
- **Instruct the patient to consult a physician** immediately if there are any signs of tissue edema and/or inflammation/infection.
- **Instruct the patient to consult a physician** if leakage through or around the voice prosthesis occurs. Leakage may cause aspiration pneumonia.
- **If used, choose laryngectomy tubes or stoma buttons** with a suitable shape that do not exert pressure on the prosthesis or catch onto the tracheal flange of the prosthesis during insertion and removal of the laryngectomy tube or stoma button. This may lead to severe tissue damage and/or accidental ingestion of the prosthesis.
- **Instruct the patient to use only genuine Provox accessories** of corresponding size (Brush, Flush, Plug) for maintenance and to avoid all other kinds of manipulation.
- **Re-use and re-processing** may cause cross-contamination and damage to the device, which could cause patient harm.

**1.6 PRECAUTIONS.**

- Carefully assess any patient with bleeding disorders or who is undergoing anticoagulant treatment for the risk of bleeding or hemorrhage prior to placement or replacement of the prosthesis

- Inspect the package before use. If the package is damaged or opened, do not use the product.
- Always use aseptic technique when handling the prosthesis in order to reduce infection risk.

## 2. Instructions for use

### 2.1 Choose size of the Voice prosthesis

#### Choosing the right shaft diameter and length of the replacement prosthesis

Be sure to use a Provox Vega voice prosthesis of the proper shaft diameter and length. Provox Vega is available in different shaft diameters and in several lengths.

- **Selecting shaft diameter**  
The clinician should determine the proper diameter of the prosthesis appropriate for the patient.  
If the selected diameter is larger than the previous prosthesis the tract must be dilated, using the Provox Dilator, appropriate to the diameter of the prosthesis being inserted.  
If a prosthesis with a smaller shaft diameter is inserted, observe and ensure that the puncture shrinks to the appropriate diameter.
- **Selecting shaft length**  
To select the correct length, you may use the current prosthesis as its own measuring device.  
If there is too much (i.e. 3 mm /~0.12 inches, or more) space between the Tracheal Flange of the old prosthesis and the mucosal wall, a shorter prosthesis should be used. If the prosthesis sits too tight, a longer prosthesis should be used.

### 2.2 Preparation

(Fig. 3-6)

#### Position voice prosthesis

1. Ensure the Voice Prosthesis is properly positioned on the Insertion Pin, firmly attached, and with the tip of the pin positioned all the way into the blue ring of the voice prosthesis (Fig. 3).

#### Fold the esophageal flange

2. Verify that the Insertion Pin is correctly positioned with the Folding Tool (The Pin shall be snapped into the Folding Tool).
3. Squeeze the Folding Tool together with two fingers (Fig. 4)
4. Attach the Loading Tube while keeping the Folding Tool closed and twist the Loading Tube until it locks in place (Fig. 5).

#### Load

5. Push the Insertion Pin forward until the voice prosthesis is aligned with the visible ring on the Insertion tube (Fig. 6).

#### Remove the old voice prosthesis

6. Remove the current (old) prosthesis from the TE-puncture by pulling it out with a non-toothed hemostat. Alternatively, at the clinician's discretion, the tracheal flange of the prosthesis can be grasped with forceps and cut off. The rest of the prosthesis is then pushed into the esophagus for passage through the intestinal tract. The patient's history of any intestinal diseases should be taken into account before using this method.

#### Prepare the puncture (optional)

7. The puncture may be dilated to prepare for the insertion of the voice prosthesis. This is usually not necessary but may facilitate insertion in patients with angled or tight punctures that easily collapse.

### 2.3 Insertion, Anterograde replacement procedure

Three methods can be used to insert a Provox Vega voice prosthesis with the Provox Insertion System. The choice of method is at the clinician's discretion.

2.3.1 System Insertion: The entire insertion system is used.

2.3.2 Tube Insertion: The folding Tool is removed before insertion.

2.3.3 Overshoot Insertion: The entire voice prosthesis is deployed in the esophagus and then retracted to the intended position.

#### 2.3.1 Method 1: System Insertion

1. **Enter the TE-puncture**  
Hold the Provox Insertion System by the Loading Tube. Enter the puncture with the Loading Tube Tip. Proceed with care if you encounter resistance. If there is resistance, dilatation and/or lubrication can ease the insertion.
2. **Insert the voice prosthesis**  
Hold the Loading Tube stable with one hand and push the Insertion Pin with the other hand until the arrows at the Distal grip surface are flush with the bottom edge of the Folding Tool (Fig. 7). At this point, the esophageal flange of the voice prosthesis is entirely unfolded in the esophagus.
3. **Release the voice prosthesis**  
Pull the Loading Tube and Folding Tool together straight out from the puncture. The voice prosthesis remains in the puncture; still firmly attached to the Insertion Pin.
4. **Finalize the procedure**  
After insertion, finalize the procedure as described below in section 2.5.

### 2.3.2 Method 2: Tube Insertion

To enhance visibility, the Folding Tool can be removed once the prosthesis has been pushed into the Loading Tube.

**Note:** It is important to always use the Provox Insertion System with Folding Tool attached when loading the prosthesis to ensure proper folding of the esophageal flange.

#### 1. Remove the Folding Tool

Remove the Folding Tool by unlocking and disconnecting it from the Loading Tube (Fig. 8).

#### 2. Enter the TE-puncture

Hold the Provox Insertion System by the Loading Tube. Enter the puncture with the Loading Tube Tip. Proceed with care if you encounter resistance. If there is resistance, dilatation and/or lubrication can ease the insertion.

#### 3. Insert the voice prosthesis

Push the Insertion Pin with the other hand until the arrows at the Proximal grip surface are flush with the Loading Tube (Fig. 9). At this point, the esophageal flange of the voice prosthesis is entirely unfolded in the esophagus.

#### 4. Release the voice prosthesis

Hold the Loading Tube with one hand and the Insertion Pin with the other hand. Pull both the Loading Tube and the Insertion Pin back until you feel the esophageal flange has reached the anterior esophageal wall. Then, let go of the Insertion Pin and pull the Loading Tube straight out from the puncture. This will release the tracheal flange.

**Note:** If the tracheal flange unfolds inside the esophagus or puncture (unintentional overshooting), grasp the tracheal flange with a non-toothed hemostat and pull/rotate the voice prosthesis into place.

#### 5. Finalize the procedure

After insertion, finalize the procedure as described below in section 2.5

### 2.3.3 Method 3: Overshoot Insertion

Overshooting of the Provox Vega voice prosthesis is occasionally performed to ensure complete development of the esophageal flange in difficult TE punctures.

#### 1. Remove the Folding Tool (Optional)

After the prosthesis has been pushed into the Loading Tube, remove the Folding Tool by unlocking and disconnecting it from the Loading Tube (Fig. 8).

#### 2. Enter the TE-puncture

Hold the Provox Insertion System by the Loading Tube. Enter the puncture with the Loading Tube Tip. Proceed with care if you encounter resistance. If there is resistance, dilatation and/or lubrication can ease the insertion.

#### 3. Insert the voice prosthesis

Hold the Loading Tube stable with one hand and push the Insertion Pin with the other hand until beyond the Proximal grip surface (Fig. 10). At this point the voice prosthesis is fully deployed in the esophagus (Fig.11).

If overshoot insertion is performed with the Folding Tool kept assembled, push the Insertion Pin beyond the Distal grip surface. At this point the voice prosthesis is fully deployed in the esophagus (Fig. 12).

#### 4. Release the voice prosthesis

Pull the Loading Tube straight out from the puncture. The prosthesis remains in the esophagus; still firmly attached to the Insertion Pin. Grasp the tracheal flange with a non-toothed hemostat and pull/rotate the voice prosthesis into place.

#### 5. Finalize the procedure

After insertion, finalize the procedure as described below in section 2.5

## 2.4 Assembly and reloading the Provox Insertion System

If the first attempt to insert the Provox Vega into the TE-puncture was unsuccessful, it is possible to reload the voice prosthesis into the Provox Insertion System.

#### Position voice prosthesis

1. Mount the voice prosthesis with the tracheal side facing down on top of the Insertion Pin.
2. Attach the Safety Strap by leading it through the Attachment Slot from either side.
3. Ensure the Voice Prosthesis is properly positioned on the Insertion Pin, firmly attached and with the tip of the pin positioned all the way into the blue ring of the voice prosthesis (Fig. 3).

#### Connect the Folding Tool

4. Guide the Insertion Pin through the opening in the Folding Tool until it clicks in place. Pull the pin until the pin is snapped into the Folding Tool (Fig. 16). The Provox Insertion System is now ready for preparation. Continue the preparation as described above in section 2.2.

**CAUTION:** Do not reload more than twice. If the Provox Vega voice prosthesis shows any signs of damage, do not use the voice prosthesis.

## 2.5 Finalize the procedure

#### After insertion:

##### Ensure proper fit

1. Ensure the proper position of the voice prosthesis by rotating and gently pulling the Insertion Pin (Fig. 17). The prosthesis is in an ideal position when the Safety Strap is pointing downwards.

##### Test proper function

2. Brush the Provox Vega with the Provox Brush. Test proper function of the prosthesis by asking the patient to speak and by observing that the prosthesis does not leak while the patient drinks water (Fig. 18).

### **Cut the Safety Strap**

3. After ensuring the voice prosthesis is placed correctly, cut off the Safety Strap so that it is flush with the Tracheal Flange (Fig. 19). The Provox Vega voice prosthesis is now ready for use.

## **3. Important patient information**

### **3.1 General information**

#### **Ensure that the patient understands to contact their clinician if:**

- There is leakage through or around the prosthesis (coughing and/or change of mucus color).
- Speaking becomes difficult (higher effort and/or voice sounds more strained).
- There are any signs of inflammation or tissue changes at the puncture site or the stoma region (pain, redness, heat, swelling, traces of blood on the brush after brushing).

#### **Also inform the patient that:**

- After a prosthesis replacement traces of blood may be found in the sputum. This may come from granulation tissue on the edges of the TE puncture.
- The Provox Vega Plug can be used to temporarily stop leakage through the voice prosthesis.
- Speech Training Sessions with a trained clinician are advised in order to acquire optimal voice sound, fluent speech, and optimal intelligibility.

### **3.2 Prosthesis maintenance**

Instruct the patient when and how to clean his or her Provox Vega voice prosthesis. Cleaning the prosthesis at least twice a day can help prolong the device life.

#### **CAUTION: Only use genuine Provox accessories that are intended for use with Provox Vega when cleaning the prosthesis.**

- Brushing the prosthesis with the Provox Brush twice a day will help to remove mucus and food remnants from the prosthesis.
- Flushing the prosthesis with the Provox Flush also helps to clean debris and mucus from the prosthesis, which can help increase the life of the device.  
**Note:** The Provox Flush is intended for use only by patients who, as assessed by the clinician who prescribes the device, have demonstrated ability to understand and consistently follow the instructions for use without clinician supervision.
- Some dietary measures, like the daily intake of yogurt or butter milk containing lactobacilli, are considered to be helpful against excessive Candida growth.

For detailed information on how to clean the Provox Vega, please see Instructions for each accessory.

## **4. Additional information**

### **4.1 Compatibility with MRI, X-ray and radiation therapy**

Provox voice prostheses (except Provox ActiValve) have been tested and found to be compatible, with Magnetic Resonance Imaging (tested up to 3 Tesla), X-ray and radiation therapy (tested up to 70 Gy). The prosthesis can be left within the TE-puncture during the examination/therapy session. Note that the radio-opaque valve seat is visible on X-ray.

### **4.2 Device lifetime**

Depending on individual biological circumstances the device life varies, and it is not possible to predict the integrity of the device over a longer period of time. The material of the device will be affected by e.g. bacteria and yeast, and the structural integrity of the device will eventually deteriorate.

Laboratory testings of simulated usage for a test period of 12 months show that, in the absence of bacteria and yeasts, the device will maintain its structural integrity for this time period. The device has not been tested for usage beyond 12 months. Usage beyond this limit is under the sole discretion of the prescriber.

### **4.3 Disposal**

Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

### **4.4 Compatibility with antifungal medicine**

Laboratory tests show only minor negative influence on the function of the Provox Vega voice prosthesis and components when using the following antifungal medications: Nystatin, Fluconazole, Miconazol.

## 5. Product model

Model	Specification
8270-18	Provox Vega 17Fr 4mm
8271-18	Provox Vega 17Fr 6mm
8272-18	Provox Vega 17Fr 8mm
8273-18	Provox Vega 17Fr 10mm
8274-18	Provox Vega 17Fr 12.5mm
8275-18	Provox Vega 17Fr 15mm
8276-18	Provox Vega 20Fr 4mm
8277-18	Provox Vega 20Fr 6mm
8278-18	Provox Vega 20Fr 8mm
8279-18	Provox Vega 20Fr 10mm
8280-18	Provox Vega 20Fr 12.5mm
8281-18	Provox Vega 20Fr 15mm
8282-18	Provox Vega 22.5Fr 4mm
8283-18	Provox Vega 22.5Fr 6mm
8284-18	Provox Vega 22.5Fr 8mm
8285-18	Provox Vega 22.5Fr 10mm
8286-18	Provox Vega 22.5Fr 12.5mm
8287-18	Provox Vega 22.5Fr 15mm

## 6. Storage conditions

Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.

## 7. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

## 8. Legal Agent and After-sales service information

Coloplast (China) Medical Devices Ltd.  
Address: Unit1301-1306, 13th Floor, Building 1, No.5  
Lido Huayuan Road, Chaoyang District,  
Beijing  
Agent contact information:  
Tel: 010-5920 1888  
Fax: 010-5920 1898  
Coloplast Customer Service Hotline  
Hotline: 400 700 7668  
Website: www.coloplast.com.cn

## 9. Production date and validity

Please see the label for the production date; the shelf life of product is 5 years.

## 10. Registration certificate number of medical devices

Certificate number of medical devices: 国械注进20223130609

## 11. Technical requirements for medical devices

Technical requirement number for medical devices: 国械注进20223130609

## 12. Edition of Instruction for Use

Article number: 11610, Edition date: 2026-03-19

## 13. Registration and Manufacturer company

Atos Medical AB  
Registration & Manufacturer location: Kraftgatan 8, SE-242 35 Hörby, Sweden  
Registration & Manufacturer contact: +46 (0)415 198 00 • info@atosmedical.com  
Production location: Kraftgatan 8, SE-242 35 Hörby, Sweden  
Country of origin: Sweden

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2026 Coloplast A/S. All rights reserved.