
PROVOX XtraHME

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
Figure 1 / 图1


Figure 2 / 图2


Figure 3 / 图3


Figure 4 / 图4

Disclaimer

Atos Medical offers no warranty—neither expressed nor implied—to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Due to local Chinese requirements, the text in the Intended use paragraph and the Overall description and product composition paragraph are not translated verbatim.

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1. Descriptive information

1.1 Intended use

The Provox XtraHME Cassette is a single use, specialized device intended for patients breathing through a tracheostoma. It is a heat and moisture exchanger (HME) that heats and humidifies inhaled air by retaining heat and moisture from exhaled air in the device. It partially restores lost breathing resistance.

1.2 Overall description and product composition

Provox XtraMoist HME and Provox XtraFlow HME are composed of a housing (Polypropylene (PP) and White PP masterbatch), Lid (PP and Beige Polyethylene (LLDPE) masterbatch) and media (Polyurethane foam polyester base with CaCl₂). Provox XtraHME Cassettes are single use devices.

1.3 CONTRAINDICATIONS

The devices should only be used in accordance with the Instructions for Use. Patients without the physical, cognitive, or mental ability required to attach, remove or operate the devices themselves, should not use the devices independently and should only use them if they are under sufficient supervision of a clinician or a trained caregiver.

Do not use the Provox XtraHME Cassette in combination with cuffed tracheal cannula; breathing may be restricted and suffocation may occur.

The devices should not be used by patients with a low tidal volume, as the added dead space may cause CO₂ (Carbon dioxide) retention, see 1.5 Technical data.

1.4 Description of the device

The Provox XtraHME Cassettes are single use devices for pulmonary rehabilitation. They are a part of the Provox HME System, which consists of HME cassettes, attachment devices and accessories.

The Provox XtraHME Cassettes have a calcium chloride treated foam in a plastic housing.

The Provox XtraHME Cassettes are available in two versions:

- Provox XtraMoist HME is intended for use during normal everyday activity.

- Provox XtraFlow HME is intended for use during physical activities since it has a lower breathing resistance. It can also be used in a two-step approach to get adapted to the higher breathing resistance of the Provox XtraMoist HME.

1.5 Technical data

Height	14.2 mm
Diameter	27.8 mm
Weight	1.5 g
Compressible Volume	5 ml (dead space)

	Provox XtraMoist	Provox XtraFlow
Pressure Drop at 30 l/min*	0.7 hPa	0.4 hPa
Pressure Drop at 60 l/min*	2.4 hPa	1.3 hPa
Pressure Drop at 90 l/min*	4.8 hPa	2.9 hPa
Moisture loss at VT=1000 ml**	21.5 mg/l	24 mg/l

* Pressure drop after 1 h according to ISO 9360.

** According to ISO 9360.

It is recommended to use the Provox XtraHME continuously. When continuously using an HME the pulmonary function is likely to improve in a majority of patients, and respiratory problems, e.g. coughing and mucus production subsequently decrease.

If you have not used HMEs previously you should be aware that the device increases breathing resistance to some extent. Especially in the beginning this may be experienced as discomfort. Starting with Provox XtraFlow Cassettes may therefore be advisable.

During the first days or weeks of use the mucus production may also appear to increase due to thinning of the mucus by retained water.

1.6 WARNINGS

- Be careful not to exert pressure on the lid of the HME unintentionally. Unintentional or accidental closing of the top lid may cause difficulty in breathing.
- Always inform the patient, caretakers and others about the closing feature of the HME cassette to ensure that they understand its function. Closing the

airway to allow voicing is a well-known feature for the laryngectomized patient with a voice prosthesis. For patients without a voice prosthesis or tracheostomized patients this feature might be unknown.

1.7 PRECAUTIONS

- Always test the function of the Provox XtraHME Cassette prior to use. The top lid should immediately return to its open position after releasing the finger.
- Do not disassemble the Provox XtraHME Cassette since this will interfere with its proper function.
- Do not reuse the Provox XtraHME Cassette or attempt to rinse it with water or any other substance. This will substantially reduce the function of the HME. Additionally the risk of potential infections may increase due to bacterial colonization of the foam.
- Do not use the Provox XtraHME Cassette longer than 24 hours. The risk of potential infections may increase with the time of use due to bacterial colonization of the foam.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device. The HME will become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.

2. Instructions for use

2.1 Operating instruction

Insert the Provox XtraHME Cassette into the connector of the attachment device (Fig. 1 or 2). Breathe normally.

To speak, press the top lid of the Provox XtraHME Cassette down with a finger (Fig. 3).

Note: Always release the lid completely at inhalation to avoid increased breathing resistance.

To remove the Provox XtraHME Cassette, hold the attachment device in place with two fingers and remove the HME Cassette from the holder (Fig. 4).

2.2 Device lifetime and disposal

The HME is for single use and must be replaced at least every 24 hours, or more often if needed.

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

2.3 Accessories

Provox XtraHME Cassettes are mainly intended to be used with the other components of the Provox HME System: the Provox Adhesives, the Provox LaryTubes and the Provox LaryButtons.

The Provox XtraHME Cassette can also be used with tracheostomy tubes that have a 22 mm connector, or with a 15 mm connector together with the Provox HME Cassette Adaptor.

For more detailed information read the Instructions for Use accompanying each product.

3. Product model

Model	Specification
7272-18	Provox XtraFlow HME (20 pcs)
7291-18	Provox XtraFlow HME (30 pcs)
8229-18	Provox XtraFlow & XtraMoist HME (5+5 pcs)
7273-18	Provox XtraMoist HME (20 pcs)
7290-18	Provox XtraMoist HME (30 pcs)

4. Additional information

4.1 Compatibility with MRI Examination

MR-Safe: This device does not contain any metallic elements and has no potential for interaction with the MRI field.

5. Storage conditions

Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C

6. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

7. Registration and Manufacturer company

Atos Medical AB

Registration & Manufacturer location: Kraftgatan 8, SE-242 35 Hörby, Sweden

Registration & Manufacturer contact: +46 (0)415 198 00 • info@atosmedical.com

Production location: Kraftgatan 8, SE-242 35 Hörby, Sweden

Country of origin: Sweden

8. Legal Agent and After-sales service information

Coloplast (China) Medical Devices Ltd.

Address: Unit1301-1306, 13th Floor, Building 1, No.5

Lido Huayuan Road, Chaoyang District,

Beijing

Agent contact information:

Tel: 010-5920 1888

Fax: 010-5920 1898

Coloplast Customer Service Hotline

Hotline: 400 700 7668

Website: www.coloplast.com.cn

9. Production date and validity

Please see the label for the production date; the shelf life of product is 3 years.

10. Registration certificate number of medical devices

Certificate number of medical devices: 国械注进 20222080154

11. Technical requirements for medical devices

Technical requirement number for medical devices: 国械注进 20222080154

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Keep away from sunlight and keep dry; 怕雨, 怕晒

Storage temperature limit; 储存温度限制

Store at room temperature. Temporary deviations within the temperature range (max-min) are allowed;
室温下存放。允许温度范围内(最高-最低)的暂时偏差。

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Instructions for use; 使用说明书

Medical Device; 医疗器械

Indicates that the product is in compliance with European legislation for medical devices; 表示本产品符合欧洲医疗器械法规

Recycling guidelines; 回收指南

XXXXX, NN
YYYY-MM-DD

XXXXX, NN = Reference number, 参考编号, 版本号

YYYY-MM-DD = Date of issue; 发布日期

11614, 03
2024-03-25

Atos
atosmedical.com

Atos Medical AB
Kraftgatan 8, SE-242 35 Hörby, Sweden
Tel: +46 (0)415 198 00 • info@atosmedical.com

Atos and the Atos logo are trademarks of Coloplast A/S.
© 2024 Coloplast A/S. All rights reserved.