Bestellinformationen	REF	Rezeptpflichtig
Provox Vega XtraSeal 22.5 Fr 4 mm	4300	Rx
Provox Vega XtraSeal 22.5 Fr 6 mm	4301	Rx
Provox Vega XtraSeal 22.5 Fr 8 mm	4302	Rx
Provox Vega XtraSeal 22.5 Fr 10 mm	4303	Rx
Provox Vega XtraSeal 22.5 Fr 12.5 mm	4304	Rx
Provox Vega XtraSeal 22.5 Fr 15 mm	4305	Rx

Evidence & Experience™

Referenzen:

1. Kress P, Schäfer P, Schwerdtfeger FP. [The custom-fit voice prosthesis, for treatment of periprothetic leakage after tracheoesophageal voice restoration]. [Article in German] Laryngorhinootologie. 2006 Jul;85(7):496-500.
2. Lewin JS, Hutcheson KA, Barringer DA, Croegaert LE, Lisec A, Chambers MS. Customization of the voice prosthesis to prevent leakage from the enlarged tracheoesophageal puncture: results of a prospective trial. Laryngoscope. 2012 Aug;122(8):1767-72.
3. Hilgers FJ, Soolsma J, Ackerstaff AH, van den Brekel MW. A thin tracheal silicone washer to solve periprothetic leakage in laryngectomies: direct results and long-term clinical effects. Laryngoscope. 2008 April;118(4):640-5.
4. Kress P, Schäfer P, Schwerdtfeger FP, Röslar S. Are modern voice prostheses better? A lifetime comparison of 749 voice prostheses. Eur Arch Otorhinolaryngol. 2014 Jan;271(1):133-40.
5. Hancock KL, Lawson NR, Ward EC. Device life of the Provox Vega voice prosthesis. Eur Arch Otorhinolaryngol. 2013 Mar;270(4):1447-53.
6. Hilgers FJ, Ackerstaff AH, Jacobi I, Balm AJ, Tan IB, van den Brekel MW. Prospective clinical phase II study of two new indwelling voice prostheses (Provox Vega 22.5 and 20Fr) and a novel anterograde insertion device (Provox Smart Inserter). Laryngoscope. 2010 Jun;120(6):1135-43.
7. Lorenz KJ, Maier H. 3. [Voice rehabilitation after laryngectomy. Initial clinical experience with the Provox-Vega® voice prosthesis and the SmartInserter® system]. [Article in German] HNO. 2010 Dec;58(12):1174-83.
8. Ward EC, Hancock K, Lawson N, van As-Brooks CJ. Perceptual characteristics of tracheoesophageal speech production using the new indwelling Provox Vega voice prosthesis: a randomized controlled crossover trial. Head Neck. 2011 Jan;33(1):13-9.
9. Hancock K, Ward E, Lawson N, van As-Brooks CJ. A prospective, randomized comparative study of patient perceptions and preferences of two types of indwelling voice prostheses. Int J Lang Commun Disord. 2012 May-Jun;47(3):300-9.

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©Atos Medical AB, 201801 REF 8009DE