



# PROVOX<sup>®</sup> HME Life


# Figures

## Provox Life™ Home/Go/Energy/Protect HME


## Provox Life™ Night HME


4.


5.


6.






7.

8.

## Disclaimer

Atos Medical offers no warranty – neither expressed nor implied – to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## Patent and trademarks

Provox<sup>®</sup> is a registered trademark owned by Atos Medical AB, Sweden.

Provox<sup>®</sup> Life<sup>™</sup> is a trademark of Atos Medical AB.

For information about protective rights (e.g. patents), please refer to our webpage [www.atosmedical.com/patents](http://www.atosmedical.com/patents)

# Contents

EN - ENGLISH.....	7
DE - DEUTSCH .....	11
NL - NEDERLANDS.....	15
FR - FRANÇAIS.....	19
IT - ITALIANO .....	23
ES - ESPAÑOL .....	27
PT - PORTUGUÊS .....	31
SV - SVENSKA.....	35
DA - DANSK .....	39
NO - NORSK .....	43
FI - SUOMI .....	47
IS - ÍSLENSKA.....	51
ET - EESTI .....	55
HU - MAGYAR.....	59
PL - POLSKI.....	63
HR - HRVATSKI.....	67
JA - 日本語 .....	71
AR - العربية.....	75

# 1. Descriptive Information

## 1.1 Intended use

Provox Life™ HMEs are single-use heat and moisture exchangers for patients breathing through a tracheostoma.

## 1.2 CONTRAINDICATIONS

The device shall not be used by patients with reduced mental or physical cognitive ability. Patients who are unable to attach or remove the device themselves, or without sufficient knowledge how to use the device, or the cognitive ability to understand the risks connected to the use, should not use the device.

The device shall not be used by patients with a low tidal volume, as the added dead space may cause CO<sub>2</sub> (Carbon dioxide) retention.

## 1.3 Description of the device

Provox Life™ HMEs are single-use devices for pulmonary rehabilitation. They come in different levels of humidification, breathing resistance and filtration that makes them suitable for different situations.

The different Provox Life™ HMEs are:

**Home:** when taking it easy.

**Go:** when you are out and about.

**Energy:** when physically active.

**Protect:** when you need protection from bacteria, virus, dust and pollen.

**Night:** when sleeping.

Provox Life™ HMEs are compatible with Provox Life™ attachments.

## 1.4 Technical data

	Pressure Drop at 30 l/min*(Pa)	Pressure Drop at 60 l/min*(Pa)	Pressure Drop at 90 l/min*(Pa)	Moisture loss at VT=1000 ml* (mg/l)
Home	60	200	450	19.5
Go	30	100	230	22.5
Energy	15	50	110	23
Protect**	55	180	320	23

	Pressure Drop at 30 l/min*(Pa)	Pressure Drop at 60 l/min*(Pa)	Pressure Drop at 90 l/min*(Pa)	Moisture loss at VT=1000 ml* (mg/l)
Night	65	210	470	18.5

HME internal volume (Dead space)\* max 13 ml

\* According to ISO 9360

\*\* Bacterial Filtration Efficiency >98%, Virus Filtration Efficiency >98%. Test method adapted from ASTM F2101.

## 1.5 WARNINGS

- Do not exert unintentional physical pressure on the device. Unintentional or accidental closing or compression of the device may obstruct breathing.
- Do not squeeze, or apply excessive force during attachment or when the device is in place. Excessive force can result in the device being pushed through the adhesive coupling, and may obstruct breathing and block airways.
- Do not use non-compatible devices or non-original devices as it may cause personal injury or damage to the device. Do not dismantle or alter the device.
- Do not re-attach a device clogged with mucus, as it may lead to obstructed breathing.
- Do not attach a device that appears slippery by mucus or other fluids, as it will be hard to attach in a correct position.

## 1.6 PRECAUTIONS

- Always test the function of the device prior to use. The lid should immediately return to its open position after releasing it. (Not applicable for Provox Life™ Night HME).
- Pay attention when the device needs to be removed or attached during the night as the night-time situation could make it more difficult to orient the device in a correct position, and it may also be difficult to detect if contaminated.
- To reduce the risk of water entering the airways, use a Provox Life™ Shower when showering.
- The device shall not be used beyond 24 hours, is not reusable and shall not be rinsed in water or in any other solutions. This will substantially reduce the function of the device and can increase the risk for infections due to bacterial colonization.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the

tracheostoma while using the device. The device will become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.

## 2. Instructions for use

### 2.1 Operating instruction

#### 2.1.1 Provox Life™ Home/Go/Energy/Protect HME

##### **To attach**

Hold the device and press lightly until the device connects to the attachment device. See fig. 1.

##### **To speak**

Press the lid to allow speech. See fig. 2.

##### **To remove**

Hold the device and pull sideways until the device disconnects. See fig. 3.

#### 2.1.2 Provox Life™ Night HME

##### **To attach**

Hold the device by the side openings between your thumb and middle finger. With your index finger on the top of the device (where the moon symbol is imprinted), press lightly until the device connects to the attachment device. Carefully rotate the device in the attachment device to verify that it has connected correctly. See fig. 4-11.

##### **To speak**

Occlude the side openings simultaneously to allow speech. See fig. 12.

##### **To remove**

Hold the device by the side openings between your thumb and middle finger, allowing your index finger to rest on the top of the device (where the moon symbol is imprinted). Pull lightly sideways until the device disconnects. See fig. 4, 5, and 13.

### 2.2 Device lifetime and disposal

- HMEs are single use devices that are designed to provide pulmonary rehabilitation both day and night, and are to be replaced in line with daily activities.
- The number and type of HMEs needed during a 24 hour period will vary based on the situations encountered and the frequency of involuntary coughing.

- Typically, one HME is used during the night and about two HMEs during the day. The HME may need to be replaced more often if clogged due to mucus secretions or other contamination. HMEs are not reusable and cannot be rinsed in water or any other solution - this washes out the special salt that is essential for HME effectiveness.
- Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

### 3. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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