VITALO

- D** Gebrauchsanweisung
- GB** Instructions for use
- F** Notice d'utilisation
- E** Instrucciones de uso
- I** Istruzioni per l'uso
- NL** Gebruiksaanwijzing
- PL** Instrukcja obsługi
- GR** Οδηγίες χρήσης

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 Packaging unit

Purpose:

Inability by laryngectomy patients to breathe through mouth and nose means a significant loss of sense of smell and taste. When used as an olfactory tube, the **VITALO** can help. The extra-corporeal airway intermittently facilitates normal ventilation of the oral and nasal mucous membranes, re-activates the olfactory nerve and taste buds, thereby making smell detection and tasting again possible. At the same time the **VITALO** enables a first-aider to administer emergency respiration to a tracheotomy or laryngectomy patient. Using **VITALO** helps reduce the risk of infection to a minimum.

Description:

The **VITALO** consists of a collapsible tube, at the ends of which a silicone mask and mouth-piece are fitted. For ventilation of the oral and nasal mucous membranes the mouth-piece is positioned firmly between the lips and the silicone mask fixed securely over the tracheostomy (Illustration 1). As happens prior to the operation, respiratory air can now enter only through the nose. The air flows through the nose, pharynx and mouth-piece into the olfactory tube and from there through the tracheostomy into the lungs. When breathing out, the air is directed in a reverse sequence. The **VITALO** is designed for use by one patient only and can be re-used. The **VITALO** can also be used to administer emergency respiration to a tracheotomy or laryngectomy patient. The silicone mask ensures firm connection of the **VITALO** to the tracheostomy to guarantee a definite air supply to the patient's lungs. The tube connects mask with mouth-piece, providing a certain distance between patient and first-aider whilst the anatomically-conform-

ing mouth-piece ensures a close fit during breathing (Illustration 2). For respiration the instructions and regulations for first-aid measures should be observed.

Cleaning:

VITALO should be cleaned in lukewarm running water without adding strong detergent.

Safety Information:

Incorrect fitting of **VITALO** can result in insufficient respiration for the patient. Inappropriate cleaning can result in infection.

Product Disposal:

Disposal of the product may be carried out only in accordance with the current national regulations for waste materials.

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