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**Finalidade:**

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☐ Packaging unit

**Purpose:**

**PRIMA-RESP® / MRI** is a respiratory tracheal tube which enables a respiratory device to be attached to the patient via the tracheostomy. It is fitted with an inflatable cuff which prevents both the escape of respiratory air and aspiration. The **PRIMA-RESP®** contains a metal spiral that prevents the tracheal tube from collapsing during use. The spiral must not be used during radiotherapy or radiological and MRI examinations.

In these cases, the **PRIMA-RESP® MRI** can be used, since the reinforcing material in this version is made of plastic.

**Description:**

Closing the tube flange clamp (illustration 1); Opening the tube flange clamp (illustration 2)

Prior to insertion of the tube, the cuff should be test-inflated using a cuff-pressure measuring device. If the cuff is intact it is released again and the **PRIMA-RESP® / MRI** tube inserted by medical staff. For easy insertion of the tube into the tracheostomy it should be lubricated with a suitable solution.

The **PRIMA-RESP® / MRI** tube comes with an obturator to make the introduction of the tube into the tracheostomy easier. It is inserted into the tube before the **PRIMA-RESP® / MRI** tube is positioned in the tracheostomy. After positioning, the obturator may be removed from the tube. With the help of the tube flange the part-length of the tube to be inserted is steplessly adjustable. Opening and closing of the flange adjustment is performed single-handed. Caution is to be exercised here to prevent the hook bending, as secure fixation of the tube flange is not ensured if this happens (see illustrations).

When the cuff is positioned securely, care must be taken to ensure that the pressure does not exceed 25 mm. Hg for long-term use. For this reason the cuff pressure should be checked regularly using a cuff-pressure measuring device. Using the uni-connector, a respiratory tube may now be attached so that the patient can be ventilated. If the tube becomes blocked due to high build-up of secretions, the respiratory tube must be detached and the patient suctioned through the **PRIMA-RESP® / MRI** tube using suction catheter. The respiratory tube can then be re-attached. When re-moving the **PRIMA-RESP® / MRI** tube the cuff must first be released. The **PRIMA-RESP® / MRI** tube is for use by one patient only and should not be used for more than 29 days. After a maximum of 6 days the **PRIMA-RESP® / MRI** tube should be removed in order to inspect the trachea to prevent – or to recognise early signs of – stenosis. Prior to re-use, the **PRIMA-RESP® / MRI** tube must be cleaned as described in Chapter „Cleaning“.

**Caution:** This product is not intended for use with lasers.

**Cleaning:**

Cleaning of the **PRIMA-RESP® / MRI** tube is performed in a saline solution at a maximum temperature of 65 °C and without chemical additives. Afterwards we recommend disinfection of the tube using a suitable solution and according to its instructions. After successful cleaning and disinfection the cuff must be checked for tightness. Finally the **PRIMA-RESP® / MRI** tube is treated with a lubricant and using a soft cloth, as described in the instructions. Failure to follow this cleaning procedure invalidates the guarantee.

**Delivery Package Contents:**

**PRIMA-RESP® / MRI** tube  
Obturator  
Tracheal tube retainer strap

**Safety Information:**

**Important:** The metal-reinforced **PRIMA-RESP®** Ventilation Tube is not suitable for use during radiotherapy or radiological and MRI examinations. The plastic-reinforced **PRIMA-RESP® MRI** is to be used for these purposes.

If the packaging is damaged the product must not be used. Failure to follow the cleaning instructions or excessive cuff pressure can result in pressure necrosis formation. If the maximum recommended usage time is exceeded, inflammation can occur and candida deposits form.