

Tracoe Mini Tracheostomy Tubes for Children

This is a digital fillable PDF.

Product variants*

Personal Information

Customer	
Name and contact details	
Patient ID / No. (anonymised)	

To order, please specify: REF + Size + TL (e.g. REF 355-2.5-C26)

REF + Size + TL	Number of units

Has been ordered before / Repeat order		
Yes <input type="checkbox"/>	Order No.:	No <input type="checkbox"/>

*Please note that the production of a variant will increase the delivery time and the product cannot be exchanged.

Requests & additional information

Please send the completed attachment to order form by email to customersupport@atosmedical.com.

Date


Signature and company stamp

Product information

Please check the appropriate boxes to indicate the desired product length and tube size:

		Size							
		2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0
Length TL (mm)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	<input type="checkbox"/>		22	22					
	<input type="checkbox"/>		24	25	30				
	<input type="checkbox"/>		26	28	33				
	<input type="checkbox"/>		28	30	36	36		40	
	<input type="checkbox"/>	26	30	32	38	39		45	50
	<input type="checkbox"/>	28	32	34	40	42	38		54
	<input type="checkbox"/>			36					
	<input type="checkbox"/>	30	34	38		45	44	50	58
	<input type="checkbox"/> *1	32	36	40	44	48	50	55	62
	<input type="checkbox"/>	34	38	42	47	50	55	60	66
	<input type="checkbox"/>	37	41	45	50	52	60	65	70
	<input type="checkbox"/>	40	44	50	52	56	65	70	
	<input type="checkbox"/>		47	55	55	58	70	75	
	<input type="checkbox"/>		50	60	60	60	75		
<input type="checkbox"/>		55		65	65				

*1 = Standard length



TL: Length along the centre line from the start of the neck flange to the lower end of the cannula

Manufacturer

Tracoe Medical GmbH
Reichelsheimer Straße 1 / 3
55268 Nieder-Olm / Germany
Tel.: +49 6136 9169-0
Fax: +49 6136 9169-200
info-tracoe@atosmedical.com
www.tracoe.com

Head office

Atos Medical AB
Hyllie Boulevard 17
SE-215 32 Malmö, Sweden
Tel: +46 (0) 415 198 00
info@atosmedical.com
atosmedical.com