

Tracoe Compress

Dressings

Product Description

There are many essential products that play an important role in tracheostomy care. Our assortment of essentials includes neckbands, dressings, products for care & cleaning as well as clothing & protectors and accessories. They fulfill a wide range of user needs and preferences, are easy to use, comfortable, discreet, offer security, and are functional.

General Information	
Product line	Dressings
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm Germany
Country of origin	Germany
Handling	
Shelf life	5 years
Maximum Useful Life	Up to 3 days

Packaging	
Packaging unit	10 pcs
Material outer packaging:	Resealable Plastic bag


Material Specifications	
Material Compress	Polyurethane
Material Coating	Aluminum* Polyurethane film
Biocompatibility	ISO 10993-5:2009 and ISO 10993-10:2010
Storage	


Certification	
Product classification acc. MDR	I
CE Marking	 0 1 2 3
GMDN	63036
Quality system certified to	EN ISO 13485
FDA registration	Exempt/E652367

PRODUCT DATA SHEET



Product Name	REF	Dimensions (L x W x H)	Weight (g)	Slit	
Tracoe Compress AL*	REF 905	9.8 x 9.0 x 2.0 mm	3,6	No	
Tracoe Compress Two-ply	REF 907-A	10.0 x 9.0 x 6 mm	5.5	No	
Tracoe Compress Two-ply, Slit	REF 907-B	10.0 x 9.0 x 6 mm	5.5	Yes	
Tracoe Compress	REF 910-A	10.0 x 9.0 x 3 mm	3.8	No	
Tracoe Compress Slit	REF 910-B	10.0 x 9.0 x 3 mm	3.8	Yes	
Tracoe Compress Coated Two-ply, Slit	REF 957	10.0 x 8.2 x 6 mm	4.5	Yes	

Tracoe Compress Coated, Slit	REF 960	10.0 x 8.2 x 3 mm	2.5	Yes	
------------------------------	---------	-------------------	-----	-----	-------------------------------------------------------------------------------------

Tracoe Compress Slit Kids	REF 961	6.3 x 6.5 x 3 mm	>1.6	Yes	
---------------------------	---------	------------------	------	-----	-------------------------------------------------------------------------------------



Tracoe Medical GmbH,
Reichelsheimer Str. 1 / 3,
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com