

Optima PercuTrach[®] Suction

Tracheostomy tube with an atraumatic inserter which ensures gentle and safe positioning of the tracheostomy tube during dilatational tracheostomy or subsequent tube changing

Transparent outer cannula with integrated 15 mm ISO connector, a cylindrical low-pressure cuff and an integrated subglottic suction line, perfectly embedded into the outer cannula to avoid pressure necrosis. Made of flexible thermosensitive material which ensure good wearing properties

The soft thin-walled low-pressure cuff adapts superbly to the trachea anatomy and ensures reliable sealing, even where low filling pressure is involved.

The atraumatic inserter is perfectly adapted for this purpose and snaps securely into the outer cannula. Moreover, it features a practical handle which incorporates elevations that prevent slipping.

The inner cannula is securely connected into the 15 mm connector of the outer cannula with a perceptible snap. Integrated grip recesses facilitates changing the inner cannula.

Two inner cannulas are included:

- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control

Important product information is lasered on to the flange, visible at a glance. All materials are biocompatible and free of latex and DEHP.



Size	ID _{oc}	ID _{ic}	OD	Length		Qty.	Ref.no	NHS Supply Chain Code
				Ln	Lb			
7	7,0mm	5,4mm	10,5mm	70,0mm	79,0mm	1 set	102454	FDH728
8	8,0mm	6,4mm	11,6mm	75,5mm	87,0mm	1 set	102455	FDH729
9	9,0mm	7,4mm	12,6mm	81,0mm	93,0mm	1 set	102456	FDH730

Includes: Neckstrap, obturator, lubricant gel and funnel connector