

Optima Suction

Tracheostomy tubes made of flexible thermosensitive material which ensure good wearing properties. Transparent outer cannula with integrated 15 mm ISO connector and a cylindrical low-pressure cuff.

Secretion that accumulates above the cuff, can be removed through the integrated subglottic suction line. This suction line is embedded perfectly into the outer cannula to avoid pressure necrosis.

The inner cannula is securely connected into the 15 mm connector of the outer cannula with a perceptible snap. Integrated grip recesses facilitates changing the inner cannula.

Two inner cannulas are included:

- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control

The soft flange has lateral eyelets for fixing the inserted tube with a neck strap. Important product information is lasered on to the flange, visible at a glance.

All materials are biocompatible, latex-free and DEHP-free.



Size	ID _{oc}	ID _{ic}	OD	Length		Qty.	Ref.no	NHS Supply Chain Code
				Ln	Lb			
7	7,0mm	5,4mm	10,5mm	70,0mm	79mm	1 set	102054	FDH762
8	8,0mm	6,4mm	11,6mm	75,5mm	87mm	1 set	102055	FDH761
10	10,0mm	8,4mm	13,6mm	87,5mm	100mm	1 set	102057	FDH764

Includes: Neckstrap, obturator and funnel connector